
# STREAM<sup>TM</sup>

Dynamic  
Neuro-Thrombectomy Net

**INSTRUCTIONS FOR USE**

# STREAM™

## Dynamic Neuro-Thrombectomy Net

READ THIS PAMPHLET CAREFULLY IN ITS ENTIRETY BEFORE USING THE DEVICE!

This document is intended to guide the operator in the use of this product. It is not intended to serve as a reference for thrombectomy techniques. Refer to all warnings, cautions, contraindications, and instructions prior to use. Please save this pamphlet for further reference.

### INDICATIONS FOR USE

- Stream is indicated for use to restore blood flow in the neurovasculature by removing thrombus for the treatment of acute ischemic stroke to reduce disability in patients with a persistent, proximal anterior circulation, large vessel occlusion, and smaller core infarcts who have first received intravenous tissue plasminogen activator (IV t-PA). Endovascular therapy with the device should be started within 6 hours of symptom onset.
- The Stream is intended to restore blood flow in the neurovasculature by removing thrombus in patients experiencing ischemic stroke within 8 hours of symptom onset. Patients who are ineligible for intravenous tissue plasminogen activator (IV t-PA) or who fail IV t-PA therapy are candidates for treatment.

### CONTRAINDICATIONS

- > Patients with a known hyper-sensitivity to nickel-titanium.
- > Patients with stenosis and/or a pre-existing stent proximal to the thrombus site that may preclude safe recovery of the Stream.
- > Patients with angiographic evidence of carotid dissection.

### CAUTIONS

- > The Stream Dynamic Neuro-Thrombectomy Net should only be used by physicians trained in interventional neuroradiology and treatment of ischemic stroke.
- > Use the Stream in conjunction with fluoroscopic visualization agents throughout the entire procedure, in accordance with standard medical practice.
- > Use the Stream in conjunction with anticoagulant agents, in accordance with standard medical practice. To prevent thrombus formation and crystal formation in the contrast medium, continuously flush the area between the guide catheter and the microcatheter and between the microcatheter and the Stream or guidewire with the appropriate solutions.
- > Check all connections to make sure that no air enters the guide catheter or microcatheter during continuous flushing.
- > Keep the Stream hydrated throughout the procedure.

### WARNINGS

- > Before use, refer to the complete warnings, contraindications and instructions that follow.
- > When using the instrument, sound medical practice must always be followed by the user.
- > The Stream is a delicate instrument and must be handled carefully. Before use and when possible, during the procedure, carefully inspect device integrity. Do not use a device that shows signs of damage. Damage may prevent the device from functioning and may cause complications.
- > Use a new microcatheter for each new Stream device.

### DESCRIPTION

The Stream Dynamic Neuro-Thrombectomy Net is a sterile, non-pyrogenic and flexible device, with a braided net, that can be adjusted to vessel diameter. The Stream is inserted and deployed under fluoroscopic visualization. The Stream is available in several sizes and models to provide coverage for vessels ranging between 1.5–6 mm in diameter. The Stream is for single-patient use and is compatible with catheters and micro-catheters of 0.017" inner diameter (Stream 17), 0.021" inner diameter (Stream and Stream XL), or larger, intended for neuro interventional procedures. The Stream handle features two operating modes that control the braided net diameter: a free mode, which allows continues adjustments and an auto-lock mode, which allows adjustments in set increments.


Figure 1 - Stream parts

Table 1 - Stream models and dimension

Model	Effective Braid Diameter (A)	Maximal Braid Length (B)	Usable Length (C)	Tip to Fluorosaver Distance (D)	Recommended Vessel Diameter	Min. Micro-catheter ID
Stream 17	0.4 – 3.5 (mm)	25 (mm)	191 (cm)	150 (cm)	1.5 – 3.5 (mm)	0.017 (inch)
Stream	0.5 – 4.5 (mm)	33 (mm)	193 (cm)	150 (cm)	1.5 – 4.5 (mm)	0.021 (inch)
Stream XL	0.5 – 6 (mm)	38 (mm)	193 (cm)	150 (cm)	4 – 6 (mm)	0.021 (inch)

- > Do not rotate or torque the Stream to avoid injury or device damage.
- > Do not attempt to re-shape the Stream's tip, as it might damage the device.
- > The Stream is provided sterile and should only be used in a sterile environment
- > Visually inspect the sterile package for breaches of packaging integrity prior to use. Verify that it has not been damaged during shipping. Do not use if package is open or damaged.
- > Do not re-sterilize and/or reuse. Reprocessing and reuse increase the risk of infection and compromise device integrity and performance.
- > Do not autoclave or expose to extreme temperatures as it may damage the Stream.
- > Do not expose the Stream to any detergents.
- > Do not use beyond "Use By" date.

### COMPLICATIONS

Possible risks include, but are not limited to: air embolism, hematoma or hemorrhage; infection; distal embolization; vessel spasm, thrombosis, dissection, or perforation; emboli; acute occlusion; ischemia; false aneurysm formation; neurological deficits including stroke; death.

### PREPARATION PROCEDURE

- Administer anti-coagulation and anti-platelet medications, as per standard institutional guidelines.
- Using angiographic radiography, determine the location and size of the occlusion.
- Select the appropriate Stream model according to the sizes and dimensions in table 1.
- Using conventional catheterization techniques and fluoroscopic visualization, position all relevant catheters (i.e. guide catheter, Distal Access Catheter, etc.) near the occlusion site. Exact location should be decided according to anatomy and occlusion properties. Connect a Rotating Hemostasis Valve (RHV) to the fitting of the catheter, and then connect it to the continuous flush line.
- Connect a second RHV to the fitting of a microcatheter and connect it to the continuous flush line.
- Using conventional catheterization techniques and fluoroscopic visualization, and aided by a suitable guidewire, advance the microcatheter within the target vessel. Anatomy permitting, position the microcatheter tip distal to the thrombus.
- Remove the guidewire from the microcatheter. Inject contrast medium through the microcatheter to visualize occlusion and distal vasculature (if relevant).
- Flush the microcatheter.
- Remove the Stream from the package hoop.
- Carefully unsheath the braided net from the loading tube.
- Check that the Stream functions properly by actuating the handle and verifying corresponding diameter changes in the braided net (see figure 2).
- Verify that both modes of handle operation are effective (Free Mode and Auto Lock). Switch between modes by rotating the top part of the handle by 90° (see Figure 3).
- Gradually decrease the diameter of the braided net while carefully re-sheathing with the loading tube.

**CAUTION:** If any resistance is detected, stop re-sheathing, identify the problem and try again. If the problem persists, use a new device.

Figure 2 - Handle actuation



Figure 3 - Handle modes

#### INSTRUCTIONS FOR USE

1. Take the Stream and verify that the tip of the device is not protruding beyond the loading tube.
2. Loosen the RHV and introduce the Stream's loading tube through the port. Seat it firmly in the microcatheter's hub.
3. Tighten the RHV to prevent back flow but be careful not to damage the Stream.
4. Using a smooth and continuous motion, advance the Stream through the loading tube until the metal portion of the actuation shaft is inserted into the microcatheter. Retract the loading tube from the RHV and retighten if necessary.
5. Advance the Stream until the fluorosaver markers reach the RHV. Then, aided by fluoroscopic visualization, advance the Stream until its tip aligns with the tip of the microcatheter.

#### WARNINGS

- > Do not allow the distal tip to exit the loading tube before fully seating the loading tube within the microcatheter hub.
- > Do not advance the Stream against resistance without careful assessment of cause. Advancing the Stream against resistance may result in vessel and/or device damage. If any resistance is met during insertion or when advancing the Stream, make sure the path through which the Stream is inserted and advanced is as straight as possible. If resistance is still prominent and cause cannot be determined, withdraw the Stream completely and try again.
- 6. If needed, reposition the microcatheter so when exposing the braided net, approximately one quarter (¼) of its length will extend distally to the occluding thrombus.
- 7. While keeping the Stream in place, Unsheathe the Stream by retracting the microcatheter until it is just proximal to the Stream's proximal marker, exposing the braided net of the Stream.

**CAUTION:** Maintain device position in the vessel when moving or repositioning the microcatheter prior to withdrawal.

8. Tighten the microcatheter's RHV
9. Inflate guide catheter balloon to occlude the vessel if using a balloon guide catheter. Alternatively, activate the DAC's aspiration, or use any standard method to support the thrombectomy procedure. Adhere to the catheter manufacturer's instructions for use.
10. While actively adjusting the size of the braided net to match vessel anatomy, slowly and gently withdraw the Stream and microcatheter into the catheter as a single unit, while applying aspiration as needed.
11. Reduce the braided net's diameter when entering the guide catheter. Reducing the diameter when entering the guide catheter will facilitate thrombus retrieval and prevent damage to the Stream.
12. Continue aspirating until the Stream and microcatheter are nearly removed from the guide catheter.

**CAUTION:** If withdrawal into the catheter is difficult, withdraw the catheter, microcatheter and the Stream through the sheath, as a single unit (deflate the guide catheter balloon if relevant). Remove the sheath, if necessary.

Disconnect the catheter RHV and fully remove the Stream, microcatheter and RHV from the catheter, as a single unit.

13. If additional recanalization attempts are needed, then:
  - a. Using the same Stream device:
    - > Rinse the device in saline solution.
    - > Gently remove any thrombus or debris from the braided net.
    - > Closely inspect the Stream for damage. If undamaged – carefully re-insert the Stream into the loading tube and repeat the aforementioned steps starting with the "Preparation Procedure" section as applicable.

**CAUTION:** A Stream device can be used in up to two (2) revascularization attempts (passes). If there is any damage do not use the device and use a new Stream device. Do not perform more than three (3) revascularization attempts at a specific vessel.

- b. Using a new Stream device – repeat the steps starting with the "Preparation Procedure" section as applicable.

#### WARNINGS

- > Do not apply excessive force when adjusting the Stream diameter. Excessive force may cause oversizing, leading to vessel damage and/or other complications.
- > When the middle marker reaches the proximal marker, it indicates that the braided net has reached its maximal diameter, do not attempt to further increase the diameter.
- > If the braided net fails to respond to handle manipulations, carefully reposition the Stream and check for responsiveness. If it continues to fail to respond, carefully withdraw the Stream for inspection.
- > Do not withdraw the Stream against resistance without careful assessment of the cause. Withdrawing the Stream against resistance may result in vessel and/or device damage. If cause cannot be determined, gently re-sheath the Stream with the microcatheter, and remove the Stream and microcatheter as a single unit.
- > If the Stream cannot be withdrawn or re-sheathed, and is unresponsive to handle manipulation, a feature can be activated to disconnect the core wire, enabling collapse of the braided net and retrieval of the device. After this feature is activated, the braided net can no longer be manipulated, and the Stream should be withdrawn and retrieved without further manipulations. This feature is available in Stream and Stream XL models only. To activate this feature, repeatedly rotate the entire handle 20-25 full turns while continuously maintaining the actuation shaft under some tension. Consequently, the core wire inside the braided net will disconnect, and the withdrawal of the device will decrease the braided net diameter and enable retrieval or re-sheathing into a catheter.
- > Do not re-use the device after the detachment feature was activated.

#### STORAGE

Store in cool, dry and dark place.

#### WARRANTY DISCLAIMER

Although this product has been manufactured under carefully controlled conditions, the manufacturer has no control over the conditions under which this product is used. The manufacturer therefore disclaims all warranties, both expressed and implied with respect to the product, but not limited to, any implied warranty of merchantability or fitness for a particular purpose. The manufacturer shall not be liable to any person or entity for any medical expenses or any direct, incidental or consequential damages caused by any use, defect, failure or malfunction of the product, whether a claim for such damages is based upon warranty, contract, tort or otherwise. No person has any authority to bind the manufacturer to any representation or warranty with respect to the product. The exclusions and limitations set out above are not intended to, and should not be construed so as to contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held illegal, unenforceable or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of this Disclaimer of Warranty shall not be affected and all rights and obligations shall be construed and enforced as if this Disclaimer of Warranty did not contain the particular part or term held to be invalid.

#### MANUFACTURER

Perflow Medical Ltd.  
4 Hatzoran St. Netanya 4250604, Israel  
Tel: +972-3-6544011  
Fax: +972-3-5546211  
www.perflow.com

#### EC REP AUTHORIZED EUROPEAN REPRESENTATIVE

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MedNet EC-REP GmbH