angiodynamics

VenaCure EVLT

Procedure Kit with Spotlight OPS Sheath

Directions For Use	2
Instrucciones de uso.....	4
Mode d'emploi.....	6
Gebrauchsanweisung	8
Istruzioni per l'uso.....	10
Gebruiksaanwijzing	12
Brugsanvisning.....	14
Οδηγίες χρήσης.....	16
Instruções de Utilização	18
Bruksanvisning.....	20
取扱説明書	22
Használati utasítás.....	24
Návod k použití	26
Instrukcje użytkowania	28
Bruksanvisning.....	30
使用说明.....	32
사용 지침	34
Kullanma Yönergeleri	36
Naudojimo instrukcijos	38
Kasutusjuhend.....	40
Instrucțiuni de utilizare	42
Návod na použitie	44
Upute za uporabu.....	46
Упатство за употреба	48
Uputstvo za upotrebu	50
Показания к применению	52

2019-07

VenaCure EVLT

Procedure Kit with Spotlight OPS Sheath

Rx ONLY

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

WARNING

Contents supplied STERILE using an ethylene oxide (EO) process. Do not use if sterile barrier is damaged. If damage is found, call your sales representative. Inspect prior to use to verify that no damage has occurred during shipping.

For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

NOTE: Read this manual and the Laser Operator's manual thoroughly prior to using the VenaCure EVLT* procedure kit. Observe all warnings, precautions and cautions noted throughout these documents. Failure to do so may result in patient complications.

ATTENTION: Insure expiration dates are still valid.

DEVICE DESCRIPTION

The VenaCure EVLT One Piece Sheath (OPS) Fiber procedure kit is a kit for use during endovascular laser procedures in the superficial venous system.

VenaCure EVLT Procedure kit contents

- | | |
|--------|---|
| PACK A | <ul style="list-style-type: none">• Percutaneous entry needle• Double ended Guidewire |
| PACK B | <ul style="list-style-type: none">• Optical fiber with distance markers (<i>SiteMarks</i>™)• Introducer sheath with fiber clamp and entry point marker |

Additional equipment required but not supplied

- Diode Laser with suitable Laser Safety glasses
- Duplex Ultrasound

In addition to the above a variety of consumable and disposable items will be required such as local anesthetic, drapes, syringes and needles. These will vary according to physician preference.

INTENDED USE / INDICATIONS FOR USE

The Laser and VenaCure EVLT Procedure Kit are intended for use in the treatment of varicose veins and varicosities associated with superficial vein reflux of the Great Saphenous vein, and with veins in the lower limbs with superficial reflux.

The D15 Plus, D30 Plus, DELTA-15, DELTA-30 and VenaCure 1470 nm Lasers and VenaCure EVLT Kits are indicated for treatment of incompetent refluxing veins in the superficial venous system.

CONTRAINDICATIONS

- Patients with thrombus in the vein segment to be treated.
- Patients with an aneurysmal section in the vein segment to be treated.
- Patients with peripheral arterial disease as determined by an Ankle-Brachial Index < 0.9.
- Patients with an inability to ambulate.
- Patients with deep vein thrombosis (DVT).
- Patients who are pregnant or breast feeding.
- Patients in general poor health.
- Other contraindications may be raised by the individual physician at the time of treatment.

Extremely tortuous vein segments may require treatment by alternative techniques (phlebectomy, sclerotherapy).

ADVERSE EVENTS

Potential complications include, but are not limited to the following:

- Hematoma
- Infection
- Paresthesia
- Phlebitis
- Pulmonary Embolism
- Skin Burns
- Thrombophlebitis
- Thrombosis
- Vessel Perforation

WARNINGS

- Treatment of a vein located close to the skin surface may result in a skin burn.
- Paresthesia may occur from thermal damage to adjacent sensory nerves.
- Tissue not targeted for treatment must be protected from injury by direct and reflected laser energy with appropriate eye and protective wear for any person present in the operating room.
- Reuse of single-use devices creates a potential risk of patient or user infections. Contamination of the device may lead to injury, illness or death of the patient.
- Reprocessing may compromise the integrity of the device and / or lead to device failure.

PRECAUTION

- Prior to and during use, avoid damaging the fiber by striking, stressing or excessive bending. Do not coil the fiber tighter than a radius of 60 mm.
- The positions of the sitemarks on the VenaCure EVLT fiber have been matched to the introducer sheath provided in the EVLT procedure kit. Alternative sheaths must not be substituted.
- Prior to and during use, avoid tight bending the introducer sheath as this can cause kinks and damage. Do not coil the sheath tighter than a radius of 60 mm.
- The entry needle provided with this kit is matched to the size of the guidewire and should not be substituted.
- Do not tighten the compression clamp on sheath until fiber is in position.

Notification of the presence of phthalates is indicated on package label.

CAUTION: This product is manufactured using polyvinyl chloride (PVC) which contains one or more phthalates, including di(2-ethylhexyl) phthalate (DEHP), which is used to soften the polymer and make it suitable for use in this application. Small quantities of DEHP may leach out into the patient during clinical use. Animal studies have shown that high-level exposure to DEHP affects male fertility and reproductive development in rats. There is currently no scientific consensus which suggests that the same effect occurs in humans.

NOTE: Patients on anticoagulant therapy are still eligible for treatment. Most physicians will only treat patients with an I.N.R. (International Normalised Ratio) at the lower therapeutic limits.

NOTE: It is sometimes helpful to nick the skin with a scalpel at the puncture site either to aid insertion of the needle, or after guidewire insertion.

- VenaCure EVLT can be performed under local anesthesia in the physician's office or as an outpatient treatment.
- VenaCure EVLT must only be performed by a qualified physician who has received specialized training in the techniques described.

HOW SUPPLIED

Store in a cool, dry, dark place. Protect from UV light. Do not use if package is opened or damaged. Do not use if labeling is incomplete or illegible.

Do not open the package until ready to use the VenaCure EVLT Procedure Kit.

EQUIPMENT HANDLING REQUIREMENTS

AngioDynamics Diode Laser:

- Prepare for use in accordance with the Laser Operator Manual.

VenaCure EVLT Procedure Kit (PACKS A & B)

- Using sterile technique, remove both Packs (A & B) from the outer packing.
- Packs A & B may be placed in the sterile field if desired

VenaCure EVLT Procedure Kit (PACK A)

- Using sterile technique open the pack, remove all contents, and inspect for damage.
- DO NOT USE if any component is damaged.
- If damage is seen, contact Customer Service or your local distributor.

VenaCure EVLT Procedure Kit (PACK B)

- Do NOT open until access has been made to the vein using the contents of PACK A
- Using sterile technique open the pack, remove all contents, and inspect for damage.
- DO NOT USE if any component is damaged.
- If damage is seen, contact Customer Service or your local distributor.

VenaCure EVLT Laser Fiber

- When required, remove the sterile fiber from the package, keeping it coiled.
- Remove the securing ties and gently uncoil the fiber. Do not remove the SMA connector cover at this time.

OPERATIONAL INSTRUCTIONS

1. Prepare and drape the leg using standard practice.
2. Use sterile technique to remove and open Pack A.
3. Place the patient in reverse Trendelenburg position and use the percutaneous needle to gain access to the vein (under local anesthesia and ultrasound guidance).

NOTE: USE ULTRASOUND GUIDANCE THROUGHOUT THE PROCEDURE to ensure that devices are correctly located.

4. Advance the guidewire through the needle until it reaches the junction with the deep vein system.
5. Remove and discard the percutaneous needle.
6. Use sterile technique to remove and open Pack B.
7. Insert the VenaCure EVLT *Spotlight OPS** sheath over the guidewire until it reaches the junction with the deep vein system.

CAUTION: DO NOT USE EXCESSIVE FORCE to introduce the sheath. Excessive force may cause damage to, or separation of the sheath.

8. Remove and discard the guidewire.
9. Advance the VenaCure EVLT optical fiber through the *Spotlight OPS* sheath until the distal SiteMark on the optical fiber is coincident with the compression clamp (back of sheath).

CAUTION: DO NOT USE EXCESSIVE FORCE to introduce the optical fiber. Use ultrasound to confirm that the *Spotlight OPS* is correctly positioned if strong resistance is encountered.

10. Use ultrasound guidance in conjunction with the echogenic tip of the *Spotlight OPS* sheath to confirm that the tip of the sheath is 1-2 cm below the junction with the deep vein system.
11. While keeping the optical fiber in its location, withdraw the sheath until the compression clamp is coincident with the proximal SiteMark on the optical fiber.
12. Tighten the compression clamp.
13. Slide the Entry Point Marker (EPM) to the tissue access site to record the length of vein treated.
14. Place the patient in the Trendelenburg position and administer tumescent anesthesia along the entire length of vein being treated.
15. Attach the optical fiber to the laser.

WARNING: Laser protective eyewear must be worn by everyone in the treatment room.

16. Use trans-illumination of the red aiming beam through the skin, and ultrasound to make final confirmation that the tip of the optical fiber is more than 1 cm from the junction with the deep vein system.
17. Activate the laser and using the distance markers along the length of the *Spotlight OPS* sheath, withdraw the sheath (with optical fiber) at a rate consistent with the venous anatomy.

CAUTION: Stop the laser firing before the optical fiber is withdrawn through the access site.

18. Instruct the patient to ambulate.

WARNING: Instruct the patient to avoid hot baths and vigorous activity for 7 days following VenaCure EVLT.

19. Dress the patient's leg and schedule subsequent follow-up according to physician preference.

WARRANTY

AngioDynamics, Inc. warrants that reasonable care has been used in the design and manufacture of this instrument. This warranty is in lieu of and excludes all other representations and warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose. Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond AngioDynamics control directly affect the instrument and the results obtained from its use. AngioDynamics obligation under this warranty is limited to the repair or replacement of this instrument and AngioDynamics shall not be liable for any incidental or consequential losses, damages or expenses directly or indirectly arising from the use of this instrument. AngioDynamics neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. AngioDynamics assumes no liability with respect to instruments reused, reprocessed or resterilized, modified or altered in any way, and makes no warranties or representations, express or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instruments.

*AngioDynamics, the AngioDynamics logo, VenaCure EVLT, EVLT, *Spotlight OPS* and SiteMarks are trademarks and/or registered trademarks of AngioDynamics Inc., an affiliate or a subsidiary.

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