varian

Embozene™

COLOR-ADVANCED MICROSPHERES

Microspheres for Embolization

Directions for Use	2
Instrucciones de uso	10
Mode d'emploi	17
Gebrauchsanweisung	24
Istruzioni per l'uso	31
Gebruiksaanwijzing	38
Instruções de Utilização	45

51100811-01

2020-03

TABLE OF CONTENTS

WARNING	3
DEVICE DESCRIPTION	3
Syringe	3
User Information	3
Contents.....	3
INTENDED USE/INDICATIONS FOR USE	3
MAGNETIC RESONANCE IMAGING	4
CONTRAINDICATIONS	4
Contraindications Specific to UFE.....	4
Contraindications Specific to Embolization for BPH.....	4
WARNINGS	4
Warnings Specific to UFE	5
Warnings Specific to UFE and Pregnancy.....	5
Warnings Specific to Embolization for BPH	5
PRECAUTIONS	5
Interaction with Pharmaceuticals.....	6
ADVERSE EVENTS	6
HOW SUPPLIED	7
Handling and Storage.....	7
OPERATIONAL INSTRUCTIONS	7
Additional Materials and Equipment	7
Catheter Compatibility.....	7
Device Preparation.....	7
Use of Prefilled Syringe	8
Table A. Design Specifications for Embozene™ Color-Advanced Microspheres.....	9
WARRANTY	9

Embozene™

COLOR-ADVANCED MICROSPHERES

Microspheres for Embolization

Rx ONLY

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

WARNING

Contents supplied STERILE using a steam process. Do not use if sterile barrier is damaged. If damage is found, call your Customer Service representative.

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

DEVICE DESCRIPTION

Embozene Microspheres are spherical, tightly calibrated, biocompatible, non-resorbable, hydrogel microspheres coated with an inorganic perfluorinated polymer (Polyzene™-F). They are available in prefilled syringes in a range of sizes suitable for embolic therapy. Embozene Microspheres are color-coded by size.

Syringe

Embozene Microspheres are supplied in a 20 ml syringe prefilled with 1 ml or 2 ml of product suspended in a non-pyrogenic, sterile transport solution of physiological saline. The total volume of Embozene Microspheres including transport solution is approximately 7 ml. Prefilled syringes of Embozene Microspheres are packaged in a sterile, sealed tray with a peel-away lid. The label indicates the specific size of the microspheres contained in the syringe (see Table A).

User Information

Contents

Quantity	Material
1	Embozene Microspheres Prefilled Syringe

INTENDED USE/INDICATIONS FOR USE

Embozene Microspheres are indicated for embolization of the following conditions:

- Hypervascular tumors
- Arteriovenous malformations
- Uterine fibroids (UFE)
- Hepatocellular carcinoma
- Benign prostatic hyperplasia (BPH)
- Tumors of head, neck, torso, and skeletal system
- Bleeding and trauma
- Pre-operative reduction of bleeding other than in the central nervous system

This device is not intended for neurovascular use.

MAGNETIC RESONANCE IMAGING

Embozene™ Microspheres are MR Safe.

CONTRAINDICATIONS

Embolization procedures shall not be performed if:

- Patient is unable to tolerate vascular occlusion procedures.
- Vascular anatomy precludes correct catheter placement or embolic injection.
- Presence or likely onset of vasospasm.
- Presence of a blood coagulation disorder that would prohibit arterial punctures.
- Presence of severe atheromatous disease that would preclude correct catheter placement.
- Presence of patent extra-to-intra-cranial anastomoses or shunts from the arterial to the venous circulation.
- Presence of collateral vessel pathways which could potentially endanger non-targeted tissue during an embolization procedure.
- Presence of any vasculature where Embozene Microspheres could pass directly into the central nervous system, central circulatory system or other non-target territories.
- Patient has high-flow arteriovenous shunt with diameter greater than the selected Embozene Microspheres.
- Patient is pregnant.
- Patient has known allergies to barium sulfate, 3-aminopropyltrialkoxysilane, polyphosphazene or IV radiopaque contrast agent.

Contraindications Specific to UFE

In addition to the general embolization contraindications, uterine fibroid embolization procedures shall not be performed if:

- Presence of suspected or active pelvic inflammatory disease.
- Presence of malignancy of the pelvic region.
- Presence of endometrial neoplasia or hyperplasia.
- Presence of submucosal fibroids with more than 50% growth into the uterine cavity.
- Presence of pedunculated serosal fibroid as the dominant fibroid(s).
- Presence of fibroids with significant collateral feeding by vessels other than the uterine arteries.

Contraindications Specific to Embolization for BPH

In addition to the general embolization contraindications, embolization procedures for benign prostatic hyperplasia shall not be performed if:

- Evidence of prostatic cancer or bladder cancer
- Urethral stricture
- Prostate size < 40 grams
- Active prostatitis
- Interest in the preservation of fertility
- Patients with renal impairment
- Peak urinary flow rate > 12 ml / sec
- Larger bladder diverticulas or stones
- Neurogenic bladder
- Detrusor failure

WARNINGS

Vascular embolization is a high risk procedure. The procedure should be performed by specialized physicians trained in vascular embolization procedures:

- Care must be taken to choose larger sized Embozene Microspheres when embolizing arteriovenous malformations with large shunts to avoid passage of the microspheres into the venous and subsequently to the pulmonary circulation.
- Extreme caution should be used for any procedures above the neck, and risk benefit assessment should be performed to avoid non-target embolization complications.

- Risks of radiation from angiography and fluoroscopy used to visualize the blood vessels during embolization, which may include a radiation burn and risks to future fertility.
- Do not use Embozene™ Microspheres in conjunction with embolization devices based on organic solvents such as ethyl alcohol or dimethyl sulfoxide (DMSO) at the same embolization site.
- Do not use ionic contrast agent with this product. Ionic contrast agents could alter the microsphere characteristics resulting in microsphere deformation and procedure failure.
- Do not use heparinized saline as this could lead to microsphere agglomeration. Agglomeration may impede microsphere delivery through the catheter or result in non-target embolization.
- Should catheter obstruction occur, remove the catheter from the patient. Do not use forceful injection, guidewires or other instruments to dislodge the blockage.

Warnings Specific to UFE

- Do not use microspheres smaller than 500 μm .
- The diagnosis of uterine sarcoma could be delayed by taking a non-surgical approach (such as UFE) to treating fibroids. It is important to pay close attention to warning signs for sarcoma (e.g., rapid tumor growth, postmenopausal with new uterine enlargement, MRI findings) and to conduct a more thorough work-up of such patients prior to recommending UFE. Recurrent or continued tumor growth following UFE should be considered a potential warning sign for sarcoma and surgery should be considered.

Warnings Specific to UFE and Pregnancy

There is no long-term data on the effects of UFE on the ability to become pregnant and carry a fetus to term, and on the development of the fetus. This procedure should only be performed on women who do not intend future pregnancy. Women who become pregnant following UFE may be at increased risk for the following:

- Postpartum hemorrhage
- Preterm delivery
- Caesarean delivery
- Malpresentation
- Abnormal placentation

Devascularization of the uterine myometrium resulting from UFE may increase the risk of uterine rupture of women who subsequently become pregnant following UFE.

Warnings Specific to Embolization for BPH

- An appropriate urological work-up should be performed on all patients (e.g., urological history and appropriate testing, such as Prostate-Specific Antigen test and, when appropriate, biopsy to rule out carcinoma).
- The effect of prostatic embolization on male fertility has not been determined. Therefore, this procedure should only be performed on men who are willing to accept the risk of future infertility.
- Three-dimensional planning imaging (e.g., magnetic resonance angiography [MRA], computed tomographic angiography [CTA]) should be performed prior to embolization on patients who have had any previous invasive treatment to the prostate (e.g., surgery, ablation, etc.) or pelvic irradiation.

PRECAUTIONS

To maintain safety, the following precautions shall be considered:

- Each package of Embozene Microspheres is intended for single patient use only. Discard any unused material.
- Physicians using Embozene Microspheres should have appropriate training and experience in a related interventional procedure.
- Do not use Embozene Microspheres if the sterile barrier, the syringe or the package appears to be opened or damaged prior to use.

- Do not use Embozene™ Microspheres that have been improperly stored or mishandled.
- The physician should carefully select the size and quantity of Embozene Microspheres according to the lesion to be treated based on the physician's education and training and currently available scientific evidence.
- Physicians must decide the most appropriate time to stop the infusion of Embozene Microspheres. Typically the artery will accept fewer Embozene Microspheres as the treatment progresses. Proximal slowing or termination of flow may indicate that the vessel or the target area is occluded by Embozene Microspheres. Careful fluoroscopic monitoring is required.
- Microsphere embolization must be performed slowly. The injection speed and manner must be controlled. Excessive injection rate may result in retrograde flow in the vessel leading to embolization of other non-target healthy tissue or organs.
- The color of the Embozene Microspheres may be visible through the skin if injected into superficial arteries.
- If arteriovenous anastomoses, branch vessels which lead away from the targeted embolization area, or emergent vessels not evident prior to embolization are present, it can lead to non-target embolization and cause severe complications for the patient.
- Embozene Microspheres smaller than 100 µm can migrate to distal anastomotic feeders and embolize circulation to distal tissue. For this reason, smaller particles have a greater likelihood of causing unwanted ischemic injury. This should be considered prior to starting the embolization procedure.
- Ischemia of tissue adjacent to the targeted area may result from post-embolization swelling. Therefore, special care should be taken to avoid such ischemia of non-tolerant, non-targeted tissue such as the nervous system.
- Consider upsizing Embozene Microspheres if angiographic appearance of embolization does not quickly appear during injection of the microspheres.
- If there are any symptoms of non-target embolization during injection, consider stopping the procedure to evaluate the possibility of shunting. Such symptoms may include changes in patient vital signs, such as hypoxia or central nervous system changes.

Interaction with Pharmaceuticals

There are no known chemical interactions between Embozene Microspheres and pharmaceuticals.

ADVERSE EVENTS

Potential adverse events associated with the use of Embozene Microspheres include, but may not be limited to:

- Allergic reaction
- Capillary bed saturation and tissue damage
- Cerebrovascular accident (CVA)
- Complications related to catheterization (e.g., hematoma at the site of entry, clot formation at the tip of the microcatheter and subsequent dislodgement, vasospasm, nerve injury, vessel trauma [e.g., dissection, perforation, rupture])
- Death
- Foreign body reactions (e.g., pain, rash, fever, inflammation)
- Hemorrhage
- Incomplete occlusion of vascular beds or territories may give rise to the possibility of post-procedural hemorrhage, development of alternative vascular pathways, recanalization, or recurrence of symptoms
- Infection
- Ischemia at an undesirable location
- Ischemic infarction
- Neurological deficits including cranial nerve palsies

- Post-embolization syndrome
- Pulmonary embolization
- Thrombosis
- Undesirable reflux, passage/migration or placement of Embozene™ Microspheres, resulting in non-target embolization
- Vessel or lesion rupture

Potential adverse events specific to UFE:

- Hysterectomy
- Infection of the pelvic region
- Premature ovarian failure (i.e., menopause)
- Temporary or permanent stopping of menstrual bleeding
- Tissue passage, fibroid sloughing, or fibroid expulsion post-UFE
- Uterine rupture
- Uterine/Ovarian necrosis
- Vaginal discharge
- Worsening of fibroid-related symptoms or the onset of new symptoms

Potential adverse events specific to embolization for BPH:

- Acute urinary retention
- Bladder ischemia or necrosis requiring surgical intervention
- Diarrhea
- Hematuria / Hematospermia
- Rectal bleeding
- Rectal or lower colonic strictures
- Sexual dysfunction or impaired fertility
- Urethral or anal burning sensation
- Urethral stricture
- Urinary incontinence
- Urinary symptoms (e.g. dysuria, urgency, frequency)

HOW SUPPLIED

Prefilled syringes of Embozene Microspheres are packaged in a sealed tray with a peel-away lid.

Sterile

Non-pyrogenic

Do not use if package is opened or damaged.

Do not use if labeling is incomplete or illegible.

Handling and Storage

Store in a cool, dry, dark place.

Product must be used prior to expiration date on label.

OPERATIONAL INSTRUCTIONS

Additional Materials and Equipment

(1) Microsphere delivery syringe

(1) 3-way stopcock

Delivery catheter

Non-ionic contrast media

Non-heparinized saline

Catheter Compatibility

Embozene Microspheres are designed to be used with a variety of catheters and microcatheters. Select a delivery catheter of appropriate size, suitable for the dimensions of the target vessels. Embozene Microspheres can tolerate temporary compression to facilitate passage through the delivery catheter. Utilize the catheter's minimum inner diameter measurement to determine catheter-to-microsphere compatibility. You may use Table A as a reference.

Device Preparation

1. Position the catheter at the desired site and perform baseline angiography to evaluate the blood supply to the lesion.

2. Carefully select the size of Embozene Microspheres according to the size of the vessel identified and catheter used.
3. Verify that the sterile packaging was not previously compromised and the catheter is not damaged.

Use of Prefilled Syringe

1. Gently swirl the contents before opening the syringe.
2. Use only non-ionic contrast agent in accordance with the contrast agent labeling with respect to dosage.

Warning: Do not use ionic contrast agent with this product. Ionic contrast agents could alter the microsphere characteristics resulting in microsphere deformation and procedure failure.

3. Add an appropriate amount of contrast agent to the product syringe to obtain homogeneous suspension and fluoroscopic visibility. The suspension is usually obtained with a mixture of 50% contrast agent and 50% water for injection or saline for injection. Contrast agent and water for injection or saline for injection solution can be added in the same proportions to obtain a more dilute suspension.
To obtain a homogeneous suspension, rotate or gently swirl the syringe every 30 seconds to agitate the microspheres and contrast agent until a homogeneous suspension is achieved. For the smaller sizes, a homogenous suspension may be achieved in less than a minute. For the larger sizes, a few additional minutes may be required to achieve a homogeneous suspension.
4. Confirm suspension prior to delivery. If microspheres have begun to settle, gently swirl or agitate to re-suspend before delivery.
5. Purge all air from the syringe.
6. Attach the 20 ml syringe to one port of the luer-lock 3-way stopcock and a 1 ml injection syringe to another port of the stopcock. Attach a delivery catheter to the remaining port on the stopcock. Ensure the stopcock is securely attached.
7. Draw the Embozene™ Microsphere mixture slowly and gently into the injection syringe to minimize the potential of introducing air into the system.
8. Under continuous fluoroscopic control, slowly infuse Embozene Microspheres into the blood stream. Always inject under free flow conditions by fully opening the stopcock. To optimize injection through the catheter, it is recommended that the syringe remains in a horizontal position during injection.

Caution: Microsphere embolization must be performed slowly. The injection speed and manner must be controlled. Excessive injection rate may result in retrograde flow in the vessel leading to embolization of non-target healthy tissue or organs.

9. Continue infusion until the desired devascularization is achieved.
10. Once the procedural endpoint is reached, wait for 5 minutes to observe whether the microspheres redistribute themselves and re-establish flow to the target. If flow is re-established, inject an additional volume of microspheres until the final procedural endpoint is achieved.

11. At the end of the infusion, remove the catheter while maintaining gentle aspiration to avoid dislodging any residual Embozene Microspheres still inside the catheter.
12. Discard any opened Embozene Microspheres prefilled syringe units.

Table A. Design Specifications for Embozene™ Color-Advanced Microspheres

Nominal Microsphere Size	Microsphere Color	Design Specifications	Compatible Catheter size Minimum Inner Diameter in (mm)
40 µm	Black	40 µm ± 10 µm	0.008 (0.2)
75 µm	Burgundy	75 µm ± 15 µm	0.008 (0.2)
100 µm	Orange	100 µm ± 25 µm	0.008 (0.2)
250 µm	Yellow	250 µm ± 50 µm	0.009 (0.2)
400 µm	Blue	400 µm ± 50 µm	0.013 (0.33)
500 µm	Red	530 µm ± 50 µm	0.016 (0.41)
700 µm	Green	700 µm ± 50 µm	0.021 (0.53)
900 µm	Purple	900 µm ± 75 µm	0.027 (0.69)
1100 µm	Gray	1100 µm ± 75 µm	0.035 (0.89)
1300 µm	Pink	1300 µm ± 75 µm	0.038 (0.97)

WARRANTY

The Manufacturer warrants that reasonable care has been used in the design and manufacture of this instrument. **This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose.** Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond the Manufacturer's control directly affect the instrument and the results obtained from its use. The Manufacturer's obligation under this warranty is limited to the repair or replacement of this instrument and the Manufacturer shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this instrument. The Manufacturer neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. **The Manufacturer assumes no liability with respect to instruments reused, reprocessed or resterilized and makes no warranties, express or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instruments.**

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All rights reserved.
**Boston Scientific Limited
**Boston Scientific Corporation
© 2020 Boston Scientific Corporation or its affiliates.