# varian

## TANDEM™

### EMBOZENE TANDEM MICROSPHERES

Microspheres for Embolization

|                                   |           |
|-----------------------------------|-----------|
| <b>Loading Guidance</b>           | <b>2</b>  |
| <b>Guía de Carga</b>              | <b>6</b>  |
| <b>Conseils de Chargement</b>     | <b>10</b> |
| <b>Leitfaden für die Aufnahme</b> | <b>14</b> |
| <b>Guida per il Caricamento</b>   | <b>18</b> |
| <b>Vulinstructie</b>              | <b>22</b> |
| <b>Orientação de Carregamento</b> | <b>26</b> |

51100857-01

2020-03

**TABLE OF CONTENTS**

**SUPPLIES** .....3

**GENERAL INFORMATION** .....3

**STORAGE INFORMATION** .....3

    Drug loaded TANDEM™ Microspheres Extended Storage .....3

**LOADING INFORMATION** .....3

**LOADING GUIDANCE** .....4

    A Step by Step Reference for Clinicians.....4

    Step 1: Hydrate Powdered Drug with WFI.....4

        Figure 1. ....4

    Step 2: Microspheres Preparation.....4

        Figure 2. ....4

    Step 3: Drug Loading .....5

        Figure 3a.....5

        Figure 3b.....5

    Step 4: Liquid Removal.....5

        Figure 4. ....5

# TANDEM™

## EMBOZENE TANDEM MICROSPHERES

### Microspheres for Embolization

TANDEM Microspheres can load and elute drugs that may be useful in the treatment of diseases in which embolization is also effective.

This Loading Guidance will provide guidelines for the loading of:

- Doxorubicin hydrochloride (HCl) (lyophilized powder)
- Irinotecan hydrochloride (HCl) (20 mg / ml solution)

#### SUPPLIES

| Quantity | Material   |
|----------|--|
| 1        | 5 ml - 10 ml water for injection (WFI)*          |
| 2        | 5 µm sterile filter needle                       |
| 1        | Waste container (e.g. empty 100 ml solution bag) |
| 1        | 18 gauge - 20 gauge sterile needle               |

\*Not required for drugs in solution

#### GENERAL INFORMATION

- Do not use liposome formulations of the drugs.
- In order to minimize the risk of microbiological contamination, TANDEM Microspheres' drug loading should be performed under controlled aseptic conditions.
- See corresponding drug package insert for dosage, handling, storage, toxicities, adverse events, and contraindications regarding its use.

#### STORAGE INFORMATION

Please refer to the drug manufacturer's package insert for specific drug information. Once TANDEM Microspheres have been loaded with drug, it is recommended they be used immediately. If storage is required, they should be stored in the dark at 2 °C - 8 °C and be used within 48 hours.

##### Drug loaded TANDEM Microspheres Extended Storage

Please refer to specific drug storage guidance for doxorubicin HCl and irinotecan HCl in the manufacturer's drug package insert.

Physical and chemical testing demonstrates that the drug loaded TANDEM Microsphere device meets device specifications enabling storage for up to 14 days when refrigerated in the dark at 2 °C - 8 °C. However, use after 48 hours is strongly discouraged due to the increased risk of microbiological contamination.

#### LOADING INFORMATION

##### Recommended drug loading dose:

- Maximum loading dose is 50 mg doxorubicin HCl and irinotecan HCl / 1 ml Microspheres

**The maximum loading dose can be achieved through the use of the following solutions:**

Doxorubicin HCl (powder):

- For 2 ml syringe, add 5 ml WFI to powder (up to 100 mg); use resulting 5 ml solution
- For 3 ml syringe, add 7.5 ml WFI to powder (up to 150 mg); use resulting 7.5 ml solution

Irinotecan HCl (solution, 20 mg / ml):

- For 2 ml syringe, use 5 ml of solution
- For 3 ml syringe, use 7.5 ml of solution

**Recommended loading times to achieve drug loading ( $98 \pm 2\%$ ):**

- Doxorubicin HCl (powder) loading time is 60 minutes
- Irinotecan HCl (20 mg / ml) loading time is 30 minutes

**LOADING GUIDANCE****A Step by Step Reference for Clinicians****Step 1: Hydrate Powdered Drug with WFI**

Skip Step 1 if drug is in solution.

- 1.1 Hydrate powdered drug with WFI - for required amount of WFI and maximum dose. See **LOADING INFORMATION** section above.


Figure 1.

**Step 2: Microspheres Preparation**

- 2.1 Open box and remove inner packaging. Verify sterile package is free of damage. Do not use if sterile barrier is damaged.
- 2.2 Peel lid to open sterile package and remove syringe from tray.
- 2.3 Gently swirl contents and allow microspheres to settle towards the plunger. It may take several minutes to fully settle.
- 2.4 Attach a sterile filter needle to the syringe. Eject transport solution.

---

**Note:** Avoid disturbing microspheres. Never compress microspheres.

---

- 2.5 Remove the filter needle and attach a sterile needle.


Figure 2.

### Step 3: Drug Loading

- 3.1 Transfer drug from its container to the TANDEM™ Microspheres syringe by drawing solution into syringe via the attached sterile needle. Leave an air bubble in the syringe for better mixing.
- 3.2 Re-cap syringe.
- 3.3 Mix by gently inverting / swirling the syringe every 5 minutes. Allow for sufficient loading time. See **LOADING INFORMATION** section above.
- 3.4 After drug loading is complete, if not used immediately, store in the dark at 2 °C - 8 °C. See **STORAGE INFORMATION** section above.


Figure 3a.


Figure 3b.

### Step 4: Liquid Removal

- 4.1 Prior to use, allow the drug loaded microspheres to settle towards the plunger.
- 4.2 Attach a sterile filter needle to the syringe and eject supernatant solution.

**Note:** Avoid disturbing microspheres. Never compress microspheres.

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All rights reserved.
**Boston Scientific Limited
**Boston Scientific Corporation
© 2020 Boston Scientific Corporation or its affiliates.