Abrax™

EN
DE
FR
ES

Sirolimus-Eluting Biodegradable Polymer Coronary Stent System - 0,014"
Biologisch abbaubares sirolimusbeschichtetes koronares Stentsystem - 0,014"
Endoprothèse coronaire en polymère biodégradable à élution de sirolimus -0,014"
Sistema de stent coronario de liberación de sirolimus con polímero biodegradable -0,014"

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ENGLISH

1. DEVICE NAME

The Abrax™ Sirolimus-Eluting Biodegradable Polymer Coronary Stent System (Abrax Stent for short) is a premounted balloon expandable drug-eluting stent with the following attributes:

A stainless steel stent with Carbon Ion Implantation on surface with a conformal coating of a biodegradable polymer carrier loaded with 2.0 µg/mm² Sirolimus in a specified release formulation with a maximum nominal content of 450 µg on the largest stent (4.0 mm x 38 mm). The proportion of Sirolimus in the coating amounts to 33% by weight. Rapid exchange balloon catheter. Two radiopaque markers which aid in the accurate placement of the stent. A balloon enabling high pressure inflations that can be used for post stent dilatation. Stent diameters of 2.0 to 4.0 mm, and stent lengths of 10 to 38 mm.

2. INDICATIONS AND USAGE

The Abrax™ Sirolimus-Eluting Biodegradable Polymer Coronary Stent System is indicated for improving coronary luminal diameter in patients with symptomatic ischemic disease due to discrete de novo and in-stent restenotic lesions (length 38mm) in native coronary arteries with vessel diameter of 2.0 mm to 4.0mm.

The Abrax™ Sirolimus-Eluting Biodegradable Polymer Coronary Stent System has been shown to significantly reduce binary restenosis, target lesion revascularisation and angiographic late loss after six and nine months when compared to Paclitaxel eluting stents.

The Abrax Stent is also indicated for treatment of abrupt or threatening closure in patients with failed interventional therapy. The treated lesion (>50%) length should be less than the nominal stent length (10 to 38 mm) with reference vessel diameters from 2.0 to 4.0 mm.

3. CONTRAINDICATIONS

The Abrax™ Sirolimus-Eluting Biodegradable Polymer Coronary Stent System is contraindicated for patients with:

- Known sensitivity to Sirolimus
- Known allergy to stainless steel
- Known allergy to PLGA polymer
- Severe reaction to contrast agents
- Patients in whom anti-platelet and/or anticoagulant therapy is contraindicated.
- In-stent Restenosis
- Myocardial infarction < 72 hours.
- Stenting of Saphenous Vein Grafts

- Unprotected left main coronary artery
- Total occlusion of target vessel
- Heavily calcified lesions
- Lesions involving arterial segments with highly tortuous anatomy
- Lesions involving a bifurcation
- Left ventricular ejection fraction < 30 % .
- Cardiogenic shock
- Presence of definite or probable intraluminal thrombus.
- Any patients judged to have a lesion which may prevent proper stent deployment

4. WARNINGS and PRECAUTIONS

Warnings

- The device carries an associated risk of subacute thrombosis, vascular complications (Complications can include bleeding, hematoma or pseudoaneurysm), and/or bleeding events. Therefore, patients should be carefully selected, and Long-life antiplatelet therapy must be prescribed.
- The device should not be used in patients with a known hypersensitivity to Sirolimus or stainless steel or poly (d-lactid-coglycolid).
- Implantation of the stent should be performed only by physicians who have received appropriate training.
- Stent placement should only be performed at hospitals where emergency coronary artery bypass graft surgery can be readily performed.
- The device should be manipulated while under high-quality fluoroscopic observation. Do not advance or retract the catheter unless the balloon is fully deflated under vacuum. If resistance is met during manipulation, determine the cause of resistance before proceeding.
- Never try to straighten a kinked hypotube. Straightening a kinked metal shaft may result in breakage of the shaft.
- Balloon pressure should not exceed the rated burst pressure. The rated burst pressure is based on the results of in vitro testing. At least 99.9% of the balloons (with 95% confidence) will not burst at or below their rated burst pressure. Use of a pressure-monitoring device is recommended to prevent over pressurization.
- Use the device before the "Expiry" date specified on the package.
- Direct Stenting without prior dilation of the lesion has not been studied using this product (note Precaution regarding predilation prior to stent implantation).
- Subsequent restenosis may require repeated dilation of the arterial segment containing the stent. The long-term outcome following repeated dilation of coronary stents is unknown at present.
- When multiple stents are required, if placement results in stent to

stent contact, stent materials should be of similar composition to avoid the possibility of dissimilar metal corrosion.

- The use of Sirolimus-Eluting stents could cause the risk of a possible inflammatory and/or prothrombotic reaction induced by the polymer coating of the stents. The responsible physician should in each case calculate the potential risk for the patient compared to the advantages of the use of a Sirolimus-Eluting stent.
- The Abrax Stent is not recommended for use in lesions that require stent overlapping.
- Due to the known risks of the implantation of drug eluting stents in combination with the following prolonged dual oral anti-platelet therapy the user must reconsider the alternative of a bypass surgery with the accompanying known risks.

Precautions

Also see **Reuse Precaution Statement**

4.1 Stent Handling – Precautions (also see Section 7 Operator's Instruction)

- Note product "Expire" date.
- Verify visually that the stent is located between the proximal and distal balloon markers.
- The Abrax™ Sirolimus-Eluting Biodegradable Polymer Coronary Stent System is designed for use as a unit. The stent is not to be removed from its delivery balloon. The stent is not designed to be crimped onto another balloon.
- Removing the stent from its delivery balloon may damage the stent and/or lead to stent embolization.
- Special care must be taken not to handle or in any way disrupt the stent position on the delivery device. This is most important during catheter removal from packaging, placement over guidewire, and advancement through hemostasis valve adapter and guiding catheter hub.
- Excessive manipulation, e.g., rolling the mounted stent, may cause coating damage or dislodgement of the stent from the delivery balloon.
- In the event the Abrax Stent is not deployed, please return the product to Rontis and avoid handling of the stent with bare hands.
- Use only the appropriate balloon inflation media (see Section "Operator's Instructions"). Do not use air or any gas medium to inflate the balloon.
- Stent contact with any fluid prior to placement is not recommended as there is a possibility of drug release. However, if it is absolutely necessary to flush the stent with sterile/isotonic saline, contact time should be limited (1 minute maximum).

4.2 Stent placement – Precautions

- Do not prepare or pre-inflate balloon prior to stent deployment other than as directed. Use balloon purging technique described in Section "Operator's Instructions".
- Implanting a stent may lead to dissection of the vessel distal and/or proximal to the stented portion, and may cause acute closure of the vessel requiring additional intervention (e.g., CABG, further dilation, placement of additional stents, or other).
- The target lesion must be sufficiently predilated prior to stent implantation.
- Do not expand the stent if it is not properly positioned in the vessel (see Section "Stent System Removal – Precautions").

- Placement of the stent has the potential to compromise side branch patency.
- Balloon pressures should be monitored during inflation. Do not exceed rated burst pressure as indicated on product label. Use of pressures higher than specified on product label may result in a ruptured balloon and potential intimal damage and dissection. The vessel should be pre-dilated with appropriate diameter balloon having a 1:1 ratio with the vessel diameter.
- Stent retrieval methods (use of additional wires, snares and/or forceps) may result in additional trauma to the vascular site. Complications can include bleeding, hematoma or pseudoaneurysm.
- If unusual resistance is felt at any time during lesion access before stent implantation, the Stent System and the guiding catheter should be removed as a single unit (see "Stent System Removal – Precautions").
- Do not attempt to pull an unexpanded stent back into the guiding catheter while is engaged in the coronary arteries, as stent damage or stent dislodgement from the balloon may occur (see "Stent System Removal – Precautions").
- An unexpanded stent should be introduced into the coronary arteries one time only. An unexpanded stent should not be subsequently moved in and out through the distal end of the guiding catheter as stent damage or stent dislodgement from the balloon may occur.

4.3 Stent System Removal – Precautions

If unusual resistance is felt at any time during lesion access before stent implantation, the Stent System and guiding catheter should be removed as a single unit.

Do not attempt to pull an unexpanded stent back into the guiding catheter while engaged in the coronary arteries, as stent damage or stent dislodgement from the balloon may occur.

When removing the entire Stent System as a single unit:

NOTE: The following steps should be executed under direct visualization using fluoroscopy.

- Maintain guidewire placement across the lesion during the entire removal process. Carefully pull back the stent system until the proximal balloon marker of the Stent System is aligned with the distal tip of the guiding catheter.
- The Stent System and guiding catheter should be pulled back until the tip of the guiding catheter is just distal to the arterial sheath, allowing the guiding catheter to straighten. Carefully retract the Stent System into the guiding catheter and remove the Stent System and guiding catheter from the patient as a single unit while leaving the guidewire across the lesion.

Failure to follow these steps and/or applying excessive force to the Stent System can potentially result in stent damage, stent dislodgement from the balloon and/or damage to the Delivery System.

4.4 Post Implant – Precautions

- Care must be exercised when crossing a newly deployed stent with an intravascular ultrasound (IVUS) catheter, a coronary guidewire, or a balloon catheter to avoid disrupting the stent geometry or coating.
- Do not perform Magnetic Resonance Imaging (MRI) scan on patient's post-stent implantation until the stent has been completely

endothelialized (90 days) to minimize the potential for migration. The stent may cause artefacts in MRI scans due to distortion of the magnetic field.

- Prescribe a long-life antiplatelet therapy

4.5 Drug Interactions

While no specific clinical data are available, drugs, like tacrolimus, that act through the same binding protein (FKBP) may interfere with the efficacy of Sirolimus.

Drug interaction studies have not been performed. Sirolimus is metabolized by CYP3A4. Strong inhibitors of CYP3A4 (e.g. ketoconazol) might cause increased Sirolimus exposure to levels associated with systemic effects, especially if multiple stents are deployed. Systemic exposure of sirolimus should also be taken into consideration if the patient is treated concomitantly with systemic immunosuppressive therapy.

4.6 Pregnancy and Lactation

This product has not been tested in pregnant women, women who wants to be pregnant, nursing mothers or in men intending to father children, effects on the developing foetus have not been studied. While there is no contraindication, the risks and reproductive effects remain unknown.

Pregnancy

In animals, Sirolimus was embryo- and foetotoxic at clinically relevant exposures. Teratogenic effects have not been observed. The potential risk for humans is unknown. The Abrax Stent should not be used in pregnant women, unless the clinical condition of the mother requires treatment with the stent.

Lactation

It is not known whether Sirolimus is excreted in human breast milk. Sirolimus is excreted in milk of lactating rats. A risk to the suckling child cannot be excluded. A decision must be made whether to discontinue breast-feeding or to implant the stent taking into account the importance of the stent to the mother.

4.7 Antiplatelet Regimen

It is mandatory to follow a recommended drug regimen for post stent medications. Guidelines for any/all recommended medications are a matter of medical research advancements. Therefore, it is not possible to recommend a specific medication herein. It is necessary to comply with local, regional, and/or internationally accepted guidelines for interventional stenting procedures in accordance with the medical standards, as stipulated by the governing authority/society for your specific area.

The physician should use best clinical judgment in determining which of the following recommendations is suitable for each specific patient case:

- Life-long single antiplatelet therapy, usually acetylsalicylic acid (ASA), is recommended after PCI with a 150–300 mg oral loading dose of acetylsalicylic acid (ASA) followed by 75–100mg daily.
- Dual Antiplatelet Therapy (DAPT) is indicated for 6 months after stent implantation. DAPT includes the above referred dose of acetylsalicylic acid (ASA) plus a clopidogrel 300–600 mg loading dose followed by 75

mg daily. After 6 months acetylsalicylic acid (ASA) should be continued indefinitely. The use of a higher maintenance clopidogrel dose (150 mg) has been proposed in patients with high thrombotic risk (e.g. in diabetes, after recurrent myocardial infarction, after early and late stent thrombosis, for complex lesions, or in life-threatening situations should occlusion occur); however, no studies have established a short- or long-term benefit of a 150 mg daily maintenance dose.

- Shorter DAPT duration (<6 months) may be considered after implantation in patients at high bleeding risk. Patients who require premature discontinuation of antiplatelet therapy secondary to significant active bleeding should be monitored carefully for cardiac events and, once stabilized, have their antiplatelet therapy restarted as soon as possible per the discretion of their treating physicians.
- DAPT may be used for more than 6 months in patients at high ischemic risk and low bleeding risk (2014 ESC/EACTS Practice Guidelines)

It is very important that the patient is compliant with the post-procedural antiplatelet recommendations. Premature discontinuation of prescribed antiplatelet medication could result in a higher risk of thrombosis, myocardial infarction or death. Prior to PCI, if a surgical or dental procedure is anticipated that requires early discontinuation of antiplatelet therapy, the interventional cardiologist and patient should carefully consider whether a drug-eluting stent and its associated recommended antiplatelet therapy is the appropriate PCI choice. Following PCI, should a surgical or dental procedure be recommended that requires suspension of antiplatelet therapy, the risks and benefits of the procedure should be weighed against the possible risk associated with premature discontinuation of antiplatelet therapy.

5. POTENTIAL ADVERSE EVENTS

Adverse events (alphabetical order) may be associated with the implantation of a coronary stent in coronary arteries, but are not limited to the following:

- Allergic reaction
- Aneurysm
- Arrhythmias
- Death
- Dissection
- Drug reactions to antiplatelet agents/anticoagulation agents/contrast medium
- Emboli, distal (tissue, air or thrombus emboli)
- Embolization, stent
- Failure to deliver the stent to intended site
- Hemorrhage
- Hypotension/Hypertension
- Infection and pain at the insertion site
- Myocardial ischemia and /or infarction
- Occlusion
- Restenosis of stented segment (greater than 50% obstruction)
- Stroke
- Thrombosis (acute, subacute or late)
- Ventricular fibrillation
- Vessel spasm
- Stent migration
- Stent Collapse
- Stent breakage or fracture may occur during implant
- During dual oral anti-platelet therapy no surgery in other fields is possible without a high risk of complications.

There may be other potential adverse events that are unforeseen at this time.

6. HOW SUPPLIED STERILE:

This device is sterilized with ethylene oxide gas. It is intended for single use only. Nonpyrogenic. Do not use if package is opened or damaged

CONTENTS:

One (1) Abrax™ Sirolimus-Eluting Biodegradable Polymer Coronary Stent System
One (1) Instruction for Use Manual
One (1) Balloon Compliance Chart

STORAGE:

Keep in a cool, dark and dry place.

DISPOSAL INSTRUCTIONS:

After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

7. OPERATOR'S INSTRUCTIONS

7.1 Inspection Prior to Use

Carefully inspect the sterile package before opening. Do not use after the "Expire" date. If the integrity of the sterile package has been compromised prior to the product "Expire" date (e.g., damage of the package), contact Rontis for return information. Do not use if any defects are noted.

NOTE:

If at any time during use of the Premounted Stent System the stainless steel proximal shaft has been bent or kinked, do not continue to use the catheter.

7.2 Materials required (not included in Stent System package)

- Select guiding catheter(s) (inner lumen min. 1.42mm) with the appropriate configuration for the coronary artery to be treated.
- 10 or 20 ml syringe
- Rotating hemostatic valve(s) with a minimum diameter of 0.096 inch/2.39mm
- Guide wire 0.014 inch/ 0.36mm
- 60% contrast medium diluted 1:1
- Sterile physiological saline solution
- Inflation device and a three-way stopcock for balloon inflation

7.3 Preparation

Packaging Removal

1. Carefully remove the delivery system from its protective tubing for preparation of the delivery system. Do not bend or kink hypotube during removal.
2. Remove the product mandrel and stent protector by grasping the catheter just proximal to the stent (at the proximal balloon bond site), and with the other hand, grasp the stent protector on the distal end and remove gently.

NOTE: Care should be taken not to kink or bend the shaft upon application.

Balloon Preparation

1. Stent contact with any fluid is not recommended, as there is a possibility of drug release. However, if it is absolutely necessary to flush the stent with saline, contact time should be limited (1 minute maximum).
2. Prepare inflation device/syringe with diluted contrast medium.
3. Attach inflation device/syringe to stopcock; attach to inflation port. Do not bend the hypotube when connecting to inflation device/syringe.
4. With tip down, orient Stent System vertically.
5. Open stopcock to Stent System; pull negative for 15 seconds; release to neutral for contrast fill.
6. Close stopcock to Stent System; purge inflation device/syringe of all air.
7. Repeat steps 4 through 6 until all air is expelled. If bubbles persist, do not use device.
8. If a syringe was used, attach a prepared inflation device to stopcock.
9. Open stopcock to Stent System.
10. Leave on neutral.

7.4 Delivery Procedure

1. Prepare the vascular access site according to standard PTCA practice.
2. Predilate the lesion/vessel with appropriate diameter balloon having a ratio of 1:1 with the diameter of the vessel.
3. Maintain neutral pressure on inflation device attached to Stent System.
4. Backload Stent System onto proximal portion of guidewire while maintaining guidewire position across target lesion.
5. Fully open rotating hemostatic valve to allow for easy passage of the stent and prevent damage to the stent.
6. Ensure guiding catheter stability before advancing the Stent System into the coronary artery. Carefully advance the Stent System into the hub of the guiding catheter, keeping the hypotube straight.

NOTE: If unusual resistance is felt before the stent exits the guiding catheter, do not force passage. Resistance may indicate a problem and use of excessive force may result in stent damage or stent is lodgement from the balloon. Maintain guidewire placement across the lesion and remove the Stent System and guiding catheter as a single unit.

7. Advance the Stent System over the guidewire to target lesion under direct fluoroscopic visualization. Utilize the proximal and distal radiopaque balloon markers as a reference point. If the position of the stent is not optimal, it should be carefully repositioned or removed (see Section 4.3 Stent System Removal – Precautions). The inside edges of the marker bands indicate both the stent edges and balloon shoulders inflated. Expansion of the stent should not be undertaken if the stent is not properly positioned in the target lesion segment of the vessel.

NOTE: If unusual resistance is felt at any time during lesion access before stent implantation, the Stent System and the guiding catheter should be removed as a single unit (see 4.3 Stent System Removal – Precautions).

- Sufficiently tighten the rotating hemostatic valve. Stent is now ready to be deployed.

7.5 Deployment Procedure

- Inflate the Stent System expanding the stent to a minimum pressure of the nominal pressure. Higher pressures may be necessary to expand the stent to optimize stent apposition against the arterial wall. Balloon pressure must not exceed rated burst pressure.
- Maintain inflation pressure for 15-30 seconds for full expansion of the stent.
- Deflate balloon by pulling negative on inflation device until balloon is fully deflated.
- Confirm stent position and deployment using standard angiographic techniques. For optimal results, the entire stenosed arterial segment should be covered by the stent. Fluoroscopic visualization during stent expansion should be used in order to properly judge the optimum expanded stent diameter as compared to the proximal and distal coronary artery diameter(s). Optimal expansion requires that the stent be in full contact with the artery wall. All efforts should be taken to assure that the stent is not underdilated.
- If stent sizing/apposition requires optimization, readvance the Stent System balloon, or another balloon catheter of the appropriate size, to the stented area using standard angioplasty techniques.
- Inflate the balloon to the desired pressure while observing under fluoroscopy. Deflate the balloon (see Balloon Compliance Chart supplied with device).
- Reconfirm stent position and angiographic result. Repeat inflations until the desired result is achieved.

7.6 Removal Procedure

- Ensure balloon is fully deflated.
- Fully open rotating hemostatic valve
- While maintaining guidewire position and negative pressure on inflation device, withdraw Delivery System.

7.1 In Vitro Information

Refer to balloon Compliance Chart supplied with device for stent inner diameter at nominal to rated burst pressure.

8. REUSE PRECAUTION STATEMENT

Rontis's obligation under this warranty is limited to the repair or replacement of this instrument and Rontis shall not be liable for any incidental or consequential loss, damage, or expense directly or indirectly arising from the use of this instrument. Rontis neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. Rontis assumes no liability with respect to instruments reused, reprocessed or resterilized and makes no warranties, express or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instrument.

LIABILITY

The product and each component of its system have been designed, manufactured, tested and packaged with all reasonable care. The warnings contained in Rontis' instructions for use are expressly to be considered as an integral part of this provision. Rontis warrants the

product until the expiry date indicated on the same. The warranty is valid provided that the use of the product was consistent with the instructions for use. Rontis disclaims any warranty of merchantability or fitness for a particular purpose of the product. Rontis is not liable for any direct, indirect, incidental or consequential damages caused by the product. Except in the case of fraud or grave fault on the Rontis part, compensation of any damage to the buyer will not, in any event, be greater than the invoice price of the disputed products. The guarantee contained in this provision incorporates and substitutes the legal guarantees for defects and compliance, and excludes any other possible liability of Rontis, however originating, from his product supplied. These limitations of liability and warranty are not intended to contravene any mandatory provisions of law applicable. If any clause of the disclaimer is considered by a competent court to be invalid or to be in conflict with the applicable law, the remaining part of it shall not be affected and remain in full force and effect. The invalid clause shall be substituted by a valid clause which best reflects Rontis' legitimate interest in limiting its liability or warranty. No person has any authority to bind Rontis' to any warranty or liability regarding the product.

CONVERSION TABLE

Contents supplied STERILE using an ethylene oxide (EO) process. Do not use if sterile seal is damaged. If damage is found, call a Rontis Representative.

9. WARRANTY

Rontis (in short: Rontis) warrants that reasonable care has been used in the design and manufacture of this instrument. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose. Handling, storage, cleaning and sterilization of this instrument, as well as other factors relating to the patient, diagnosis, treatment, surgical procedures, and other matters beyond Rontis's control directly affect the instrument and the results obtained from its use.

CONVERSION CHART

1cc	1mL		
1French	0,0130"	0,33 mm	
1bar	0,987 atm	14,5 PSI	10 ⁵ Pa

GRAPHICAL SYMBOLS FOR MEDICAL DEVICE LABELLING

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Rontis Corporation S.A.