# Interlock™ - 35

## DETACHABLE COIL

### Fibered IDC™ Occlusion System

<b>Directions for Use</b>	<b>2</b>
<b>Instrucciones de uso</b>	<b>8</b>
<b>Mode d'emploi</b>	<b>15</b>
<b>Gebrauchsanweisung</b>	<b>22</b>
<b>Istruzioni per l'uso</b>	<b>29</b>
<b>Gebruiksaanwijzing</b>	<b>36</b>
<b>Instruções de Utilização</b>	<b>43</b>

<b>TABLE OF CONTENTS</b> .....	<b>2</b>
<b>WARNING</b> .....	<b>3</b>
<b>DEVICE DESCRIPTION</b> .....	<b>3</b>
Contents .....	<b>3</b>
Figure 1. Interlock™ - 35 Fibered IDC™ Occlusion System.....	<b>3</b>
<b>INTENDED USE/INDICATIONS FOR USE</b> .....	<b>3</b>
<b>CONTRAINDICATIONS</b> .....	<b>3</b>
<b>WARNING</b> .....	<b>3</b>
<b>PRECAUTIONS</b> .....	<b>3</b>
<b>ADVERSE EVENTS</b> .....	<b>4</b>
<b>HOW SUPPLIED</b> .....	<b>4</b>
Handling and Storage .....	<b>4</b>
<b>OPERATIONAL INSTRUCTIONS</b> .....	<b>4</b>
Figure 2. Flush device prior to use .....	<b>4</b>
Figure 3. Example of Hand Injection Setup.....	<b>5</b>
Figure 4. Example of Continuous Flush Setup .....	<b>5</b>
Figure 5. Introducer Sheath with Twist-Lock Mechanism .....	<b>5</b>
Figure 6. Unlock via Counter-Clockwise Rotation .....	<b>5</b>
Figure 7. Advancing Interlock - 35 Fibered IDC Occlusion System to Pre-Release Position.....	<b>6</b>
Figure 8. Delivering Interlock - 35 Fibered IDC Occlusion System .....	<b>6</b>
Interlock - 35 Fibered IDC Occlusion System Removal Procedure .....	<b>6</b>
Figure 9. Lock via clockwise rotation.....	<b>7</b>
<b>MAGNETIC RESONANCE IMAGING (MRI)</b> .....	<b>7</b>
3.0 Tesla Temperature Information .....	<b>7</b>
1.5 Tesla Temperature Information .....	<b>7</b>
Image Artifact Information .....	<b>7</b>
<b>WARRANTY</b> .....	<b>7</b>

# Interlock™ - 35

## DETACHABLE COIL

### Fibered IDC™ Occlusion System

#### ⚠ ONLY

**Caution:** Federal Law (USA) restricts this device to sale by or on the order of a physician.

#### WARNING

Contents supplied STERILE using an ethylene oxide (EO) process. Do not use if sterile barrier is damaged. If damage is found, call your Boston Scientific representative.

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

#### DEVICE DESCRIPTION

The Interlock - 35 Fibered IDC Occlusion System includes a coil manufactured from platinum-tungsten alloy that is mechanically attached to a coil delivery wire. This assembly is contained within an introducer sheath. The platinum coil contains synthetic fibers for greater thrombogenicity. The Interlock - 35 Fibered IDC Occlusion System is designed to be delivered under fluoroscopy through a 5F (1.70 mm) OD (0.035 in [0.89 mm] or 0.038 in [0.97 mm] inner lumen) Imager™ II Selective Diagnostic Catheter without side flushing holes. The interlocking delivery wire design allows the coil to be advanced and retracted before final placement in the vessel, thus aiding in more controlled delivery including the ability to withdraw the coil prior to deployment.

#### Contents

Quantity	Material
1	Interlock - 35 Fibered IDC Occlusion System
1	Rotating Hemostatic Valve (RHV)


Figure 1. Interlock - 35 Fibered IDC Occlusion System

#### INTENDED USE/INDICATIONS FOR USE

The Interlock - 35 Fibered IDC Occlusion System is indicated for obstructing or reducing blood flow in the peripheral vasculature during embolization procedures. This device is not intended for neurovascular use.

#### CONTRAINDICATIONS

None known.

#### WARNING

Do not advance the delivery wire once the coil has been placed. Perforation or damage of the vessel wall could occur.

#### PRECAUTIONS

Carefully read all instructions prior to use. Observe all warnings and precautions noted throughout these and other instructions relevant to the procedure. Failure to do so may result in complications.

Selection of a coil diameter smaller than the vessel diameter may result in coil migration. Coil selection is a matter of physician preference and the clinical situation however it is recommended the physician oversize the coil, as applicable based on the placed location. The shape and diameter of the vessel to be occluded as well as proximity to branch vessels generally govern selection of the coil diameter and length. During continuous flush or hand injection setup check that all fittings are secure so that air is not introduced into the catheters.

Physicians should exercise their clinical judgment in selection and use of catheters. Manufacturers may make changes to their catheters without notice which may impact their suitability for use with Interlock - 35 Fibered IDC Occlusion System. Boston Scientific Corporation provides no warranty for use of third party catheters with its products. The use of other diagnostic catheters may result in an inability to deliver, deploy, or recapture the device.

Do not attempt to use the Interlock - 35 Fibered IDC Occlusion System with a soft-walled delivery catheter, such as the Terumo Glidecath™ Catheter or AngioDynamics Soft-Vu® Catheter. The Interlock - 35 Fibered IDC Occlusion System will encounter significant resistance when advancement through a soft-walled delivery catheter is attempted.

During device preparation, do not remove the Interlock - 35 Fibered IDC Occlusion System assembly from the introducer sheath. Removal may result in pre-mature deployment of the coil.

Do not apply excessive force while seating the introducer sheath in the catheter hub. Introducer sheath tip deformation and difficulty with coil delivery into the catheter could result.

Ensure that the introducer sheath remains firmly seated in the catheter hub to prevent premature deployment during introduction of coil into catheter.

If significant friction is encountered during coil advancement, vigorously inject heparinized saline through sideport of RHV. If flushing does not rapidly resolve friction, immediately remove coil from delivery catheter to prevent damage to the coil.

Do not advance or retract the Interlock - 35 Fibered IDC Occlusion System too quickly or against significant resistance. Forcing the Interlock-35 system forward or backward against significant resistance could result in stretching of the coil or damage to the interlocking mechanism resulting in a loss of functionality. If friction is noted in any successive coil, carefully examine both coil and diagnostic catheter for possible damage. Replace both if necessary.

Do not rotate the delivery wire more than one turn (360 degrees) during delivery of the Interlock - 35 Fibered IDC Occlusion System. Excessive rotation of the delivery wire may damage the Interlock™ - 35 Fibered IDC™ Occlusion System or may result in premature detachment of the interlocking arms within the catheter.

**Caution:** Advance and retract the Interlock - 35 Fibered IDC Occlusion System smoothly, especially in tortuous anatomy. In the event resistance is felt when the coil is in the catheter, retract the system backwards until resistance decreases and then continue advancing forward again. If the Interlock-35 system will not advance forward remove and replace the entire system.

Do not advance the Interlock - 35 Fibered IDC Occlusion System if it becomes lodged within the catheter. Determine the cause of the resistance and replace the catheter and coil if necessary. See Interlock - 35 Fibered IDC Occlusion System Removal Procedure for further instructions.

Multiple embolization procedures may be required to achieve the desired occlusion of some vessels.

Replace catheters periodically during delivery of multiple coils or if increased resistance is noted during coil delivery.

Axial compression or tension forces may be stored in the catheter shaft during Interlock - 35 Fibered IDC Occlusion System delivery, and coil release may lead to catheter tip movement. During the procedure and prior to Interlock - 35 detachment, verify repeatedly that the distal shaft of the catheter is not under stress by slightly repositioning the catheter, delivery wire, or entire assembly simultaneously.

#### ADVERSE EVENTS

The complications that may result from a peripheral embolization procedure include, but are not limited to:

- Complications related to catheterization (e.g., hematoma at the site of entry, clot formation at the tip of the catheter and subsequent dislodgement, nerve and vessel dissection or perforation, etc.)
- Pain
- Hemorrhage
- Infection necessitating medical intervention
- Foreign body reactions necessitating medical intervention
- Emboli
- Ischemia
- Vasospasm
- Tissue necrosis
- Undesirable clot formation of the vasculature
- Recanalization
- Death
- Temporary neurological deficit

#### HOW SUPPLIED

The Interlock - 35 Fibered IDC Occlusion System includes an embolic coil with interlocking delivery wire, introducer sheath, and rotating hemostatic valve (RHV).

Do not use if package is opened or damaged.

Do not use if labeling is incomplete or illegible.

#### Handling and Storage

Store in a cool, dry, dark place.

#### OPERATIONAL INSTRUCTIONS

The Interlock - 35 Fibered IDC Occlusion System is recommended for use with a 5F (1.70 mm) OD (0.035 in [0.89 mm] or 0.038 in [0.97 mm] inner lumen)

Imager™ II Selective Diagnostic Catheter without side flushing holes.

In order to achieve optimal performance of the Interlock - 35 Fibered IDC Occlusion System and reduce the risk of thromboembolic complications, it is critical that the 5F Imager II Selective Diagnostic Catheter is flushed vigorously, either utilizing a continuous flush or hand injection setup, before and after the introduction of each Interlock - 35 Fibered IDC Occlusion System.

**Caution:** Physicians should exercise their clinical judgment in selection and use of catheters. Manufacturers may make changes to their catheters without notice which may impact their suitability for use with Interlock - 35 Fibered IDC Occlusion System. Boston Scientific Corporation provides no warranty for use of third party catheters with its products. The use of other diagnostic catheters may result in an inability to deliver, deploy, or recapture the device.

**Caution:** Do not attempt to use the Interlock - 35 Fibered IDC Occlusion System with a soft-walled delivery catheter, such as the Terumo Glidecath™ Catheter or AngioDynamics Soft-Vu® Catheter. The Interlock - 35 Fibered IDC Occlusion System will encounter significant resistance when advancement through a soft-walled delivery catheter is attempted.

1. Prior to use, ensure that sterile packaging is intact. Return the device if sterility appears to have been compromised.
2. Place the catheter in the area to be embolized per standard technique. Take care to position the catheter tip parallel with, not perpendicular to, the vessel wall to facilitate deployment of the coil.
3. Fill a 10cc syringe with heparinized saline solution. Attach the filled 10cc Syringe to the flush port on the dispenser hoop. Flush the dispenser hoop vigorously, bathing the coil within.



Figure 2. Flush device prior to use.

4. Slowly withdraw the Interlock - 35 Fibered IDC Occlusion System from its dispenser hoop and inspect assembly. Discard if there is any evidence of damage. Ensure that the Interlock - 35 Fibered IDC Occlusion System arms are interlocked inside the introducer sheath.


**Caution:** Do not remove the Interlock™ - 35 Fibered IDC™ Occlusion System assembly from the introducer sheath.

5. Attach the included RHV to the proximal luer adapter on the hub of the catheter. Begin setup for either (a) Hand Injection or (b) Continuous Flush of heparinized saline solution by doing the following:
  - a. Hand Injection Setup: Connect a 20cc syringe filled with heparinized saline solution to the side arm of the RHV.

**Figure 3. Example of Hand Injection Setup**

- b. Continuous Flush Setup: Attach a line for continuous flush of heparinized saline solution to the RHV. In general, one drop of saline solution every 1-3 seconds is recommended.

**Figure 4. Example of Continuous Flush Setup**


**Caution:** Check that all fittings are secure so that air is not introduced into the catheter during flushing.

**Caution:** In order to reduce the risk of complications, hand injection or continual introduction of heparinized saline solution should be maintained through the catheter and any intraluminal device. Heparinized saline solution reduces retrograde flow of blood into the catheter during coil delivery and reduces the potential for premature coil thrombosis, contrast crystal formation and/or clotting on both the coil and inside the catheter lumen.


6. Open the thumbscrew of the RHV and carefully insert the Interlock - 35 Fibered IDC Occlusion System until the distal tip of the introducer sheath is firmly seated in the catheter hub.

**Caution:** Do not apply excessive force while seating the introducer sheath in the catheter hub. Introducer sheath tip deformation and difficulty with coil delivery into the catheter could result.

7. Tighten the RHV thumbscrew to prevent retrograde flow but not so tight as to pinch the introducer sheath and inhibit forward movement of the delivery wire.
8. Prior to advancing the Interlock - 35 Fibered IDC Occlusion system from its introducer sheath into the catheter, ensure that blood is minimized within the RHV or on the fibers of the coil which are visible within the coil introducer. If blood is present, take one of the following actions:
  - a. Hand Injection Setup: Inject flush solution into RHV using the 20cc syringe connected to RHV sideport. Continue flushing until blood is minimized on coil or within RHV.
  - b. Continuous Flush Setup: Maintain in-line pressure of the continuous flush to prevent retrograde flow onto the coil. If bleed-back continues, increase the rate of continuous flush infusion until blood is minimized on the coil and within the RHV.
9. Release the Interlock - 35 Fibered IDC Occlusion System inside its introducer sheath by gently pinching the sheath on both sides of the twist-lock mechanism and rotating proximal side counter-clockwise (Figure 5 and 6).

**Figure 5. Introducer Sheath with Twist-Lock Mechanism**

**Figure 6. Unlock via Counter-Clockwise Rotation**

10. Transfer the Interlock - 35 Fibered IDC Occlusion System and delivery wire from the introducer sheath into the catheter by advancing the delivery wire in a smooth, continuous manner.


**Caution:** Ensure that the introducer sheath remains firmly seated in the catheter hub to prevent premature deployment.

**Caution:** If significant friction is encountered during coil advancement, vigorously inject heparinized saline through sideport of RHV. If flushing does not rapidly resolve friction, immediately remove coil from delivery catheter to prevent damage to the coil.

11. Gently withdraw and remove the introducer sheath once the coil is visibly transferred completely into the delivery catheter. The Interlock™ - 35 coil and interlocking arms should not be visible. Do not discard the sheath in case it is necessary to remove the Interlock - 35 Fibered IDC™ Occlusion System prior to deployment.

The following steps describe the Coil Delivery of the Interlock - 35 Fibered IDC Occlusion System.

12. Maneuver the Interlock - 35 Fibered IDC Occlusion System under fluoroscopy until the coil detachment zone is approximately 1 cm proximal to the catheter tip (Figure 7).

**Figure 7. Advancing Interlock - 35 Fibered IDC Occlusion System to Pre-Release Position**


13. If Interlock - 35 Fibered IDC Occlusion System repositioning is necessary, gently retract the Interlock - 35 Fibered IDC Occlusion System under fluoroscopy. If repositioning is difficult or impossible, remove and discard the Interlock - 35 Fibered IDC Occlusion System. (See Interlock - 35 Fibered IDC Occlusion System Removal Procedure for further instructions.)

**Caution:** Do not advance or retract the Interlock - 35 Fibered IDC Occlusion System too quickly or against significant resistance. Forcing the Interlock-35 system forward or backward against significant resistance could result in stretching of the coil or damage to the interlocking mechanism resulting in a loss of functionality. If friction is noted in any successive coil, carefully examine both coil and diagnostic catheter for possible damage. Replace both if necessary.

14. To deploy the coil, slowly advance the delivery wire under fluoroscopy until interlocking arms pass catheter's tip (Figure 8). Note: If deployment resistance is encountered, slowly rotate delivery wire until coil deploys.

**Caution:** Do not rotate the delivery wire more than one turn (360 degrees) during delivery of the Interlock - 35 Fibered IDC Occlusion System. Excessive rotation of the delivery wire may damage the Interlock - 35 Fibered IDC Occlusion System or may result in premature detachment of the interlocking arms within the catheter.

**Caution:** Advance and retract the Interlock - 35 Fibered IDC Occlusion System smoothly, especially in tortuous anatomy. In the event resistance is felt when the coil is in the catheter, retract the system backwards until resistance decreases and then continue advancing forward again. If the Interlock-35 system will not advance forward remove and replace the entire system.

**Figure 8. Delivering Interlock - 35 Fibered IDC Occlusion System**

**Warning:** Do not advance the delivery wire once the coil has been placed. Perforation or damage of the vessel wall could occur.

**Caution:** Multiple embolization procedures may be required to achieve the desired occlusion of some vessels. To position another Interlock - 35 Fibered IDC Occlusion System, return to Steps 1-15 in Directions for Use.

**Caution:** Replace catheters periodically during delivery of multiple coils or if increased resistance is noted during coil delivery.

**Caution:** Axial compression or tension forces may be stored in the catheter shaft during Interlock - 35 Fibered IDC Occlusion System delivery, and coil release may lead to catheter tip movement. During the procedure and prior to Interlock - 35 coil detachment, verify repeatedly that the distal shaft of the catheter is not under stress by slightly repositioning the catheter, delivery wire, or entire assembly simultaneously.


15. Carefully remove the delivery wire after coil deployment so that the delivery arm does not catch on the valve in the RHV thumbscrew.

#### **Interlock - 35 Fibered IDC Occlusion System Removal Procedure**

An Interlock - 35 Fibered IDC Occlusion System must be removed if the coil is determined to be the incorrect size. If resistance is noted and repositioning is difficult, remove the coil and discard. Make sure that the introducer sheath twist-lock mechanism is disengaged to facilitate threading of the sheath over the delivery wire (Figure 6).

1. Gently begin to retract the Interlock - 35 Fibered IDC Occlusion System under fluoroscopy. If resistance is encountered, retract the catheter and delivery wire simultaneously to facilitate movement.
2. Once the Interlock - 35 Fibered IDC Occlusion System is withdrawn to approximately the mid-shaft point in the catheter, gently thread the distal end of the introducer sheath over the proximal end of the delivery wire.
3. Open the thumbscrew of the RHV and carefully advance the introducer sheath until it is firmly seated in the proximal luer adapter of the catheter.

4. Tighten the RHV thumbscrew just enough to prevent retrograde flow but not so tight as to inhibit backward movement of the delivery wire through the catheter.
5. Holding the introducer sheath in place, gently withdraw the Interlock™ - 35 Fibered IDC™ Occlusion System until the interlocking arms and distal coil tip are visible inside the sheath.
6. Lock the Interlock - 35 Fibered IDC Occlusion System into position by gently pinching the introducer sheath on both sides of the twist-lock mechanism and rotating the proximal side clockwise (Figure 9).

**Figure 9. Lock via clockwise rotation**

7. Withdraw the introducer sheath/delivery wire assembly from the catheter/ RHV assembly.

### MAGNETIC RESONANCE IMAGING (MRI)

Non-clinical testing has demonstrated that the Interlock - 35 Fibered IDC Occlusion System is MR Conditional. It can be scanned safely under the following conditions:

- Static magnetic field of 1.5 Tesla or 3.0 Tesla
- Static magnetic field gradient < 25 T/m.
- Product of static magnetic field and static magnetic field gradient < 50 T<sup>2</sup>/m (extrapolated)
- Normal operating mode of the MR system with a maximum whole body averaged specific absorption rate of 2 W/kg and use of transmit/receive head coil and/or whole body transmit coils

The Interlock - 35 Fibered IDC Occlusion System should not migrate in this MRI environment. Non-clinical testing at field strengths other than 1.5 Tesla or 3 Tesla has not been performed to evaluate coil migration or heating.

### 3.0 Tesla Temperature Information

Non-clinical testing of RF-induced heating was performed at 128 MHz in a 3.0 Tesla Magnetom Trio™ Device, Siemens Medical Solutions MR system, software version Numaris/4, syngo MR A30. The coils tested were in a location and orientation in the phantom that produced the worst case Radio Frequency (RF) heating. RF power was applied for 15 minutes with the conductivity of the phantom material about 0.24 S/m. The phantom average SAR calculated using calorimetry was 3.3 W/kg. Predicted in-vivo heating based on these non-clinical tests and computer simulation of the patient exposure to the electromagnetic fields in MRI yielded the following maximal in vivo rises:

- For vessels in the body the calculated temperature rise was 4.7°C with an uncertainty upper bound temperature of 6.4°C for a whole body average SAR value of 2.0 W/kg and a continuous scan time of 15 minutes.
- The actual in vivo rise is expected to be less than these values as the calculations did not include the cooling effects due to blood flow around the coil and blood perfusion in the tissue outside the coil.

### 1.5 Tesla Temperature Information

Non-clinical testing of RF-induced heating was performed at 64 MHz in a 1.5 Tesla Intera™ Medical Magnetic Resonance Imaging Apparatus Philips Medical Systems, software version Release 10.6.2.0, 2006-03-10 whole body coil MR scanner. The coils tested were in a location and orientation in the phantom that produced the worst case RF heating. RF power was applied for 15 minutes with the conductivity of the phantom material about 0.26 S/m. The phantom average SAR calculated using calorimetry was 3.6 W/kg. Predicted in-vivo heating based on these non-clinical tests and computer simulation of the patient exposure to the electromagnetic fields in MRI yielded to the following maximal in vivo rises:

- For vessels in the body the calculated temperature rise was 5.3°C with an uncertainty upper bound temperature of 7.2°C for a whole body average SAR value of 2.0 W/kg and a continuous scan time of 15 minutes.
- The actual in vivo rise is expected to be less than these values as the calculations did not include the cooling effects due to blood flow around the coil and blood perfusion in the tissue outside the coil.

### Image Artifact Information

MR image at 1.5 and 3 Tesla may be performed immediately following the implantation of the Interlock - 35 Fibered IDC Occlusion System.

MR image quality may be compromised if the area of interest is in the same area or relatively close to the position of the coil. MR image artifact has been evaluated at 1.5 and 3 Tesla only.

The image artifact extended as much as 8 mm from the metal of the device when scanned in non-clinical testing using a Spin Echo sequence. With a Gradient Echo sequence, the image artifact extended as much as 14 mm beyond the metal of the device. Image tests were made in a 3.0 Siemens Magnetom Trio Device, software version Numaris/4. This testing was completed using ASTM F2119-07 test method.

### WARRANTY

Boston Scientific Corporation (BSC) warrants that reasonable care has been used in the design and manufacture of this instrument. **This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose.** Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond BSC's control directly affect the instrument and the results obtained from its use. BSC's obligation under this warranty is limited to the repair or replacement of this instrument and BSC shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this instrument. BSC neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. **BSC assumes no liability with respect to instruments reused, reprocessed or resterilized and makes no warranties, express or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instruments.**

Glidecath is a trademark of Terumo Corporation.

Soft-Vu is a trademark of AngioDynamics Inc.

Magnetom Trio is a trademark of Siemens.

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Intera is a trademark of Koninklijke Philips Electronics N.V.
© 2014 Boston Scientific Corporation or its affiliates. All rights reserved.
Boston Scientific Limited
Boston Scientific (Australia) Pty Ltd