Renegade™ HI-FLO™ FATHOM™ System

WITH FATHOM™ - 16 GUIDEWIRE

Pre-Loaded Microcatheter System

Directions for Use	2
Instrucciones de uso	6
Mode d'emploi	10
Gebrauchsanweisung	14
Istruzioni per l'uso	18
Gebruiksaanwijzing	22
Instruções de Utilização	26

TABLE OF CONTENTS

WARNING	3
DEVICE DESCRIPTION	3
COMPATIBILITY INFORMATION	3
Table 1. Renegade™ HI-FLO™ Compatibilities	3
Table 2. FATHOM™-16 Compatibilities	3
Contents	3
INTENDED USE/INDICATIONS FOR USE	3
CONTRAINDICATIONS	3
ADVERSE EVENTS	3
WARNINGS	4
PRECAUTIONS	4
Figure 1. Example of Continuous Flush Setup	4
PREPARATIONS FOR USE	4
DIRECTIONS FOR USE	5
HOW SUPPLIED	5
Handling and Storage	5
WARRANTY	5

Renegade™ HI-FLO™ FATHOM™ System

WITH FATHOM™ - 16 GUIDEWIRE

Pre-Loaded Microcatheter System

Rx ONLY

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

WARNING

Contents supplied STERILE using an ethylene oxide (EO) process. Do not use if sterile barrier is damaged. If damage is found, call your Boston Scientific representative.

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

DEVICE DESCRIPTION

The Renegade HI-FLO FATHOM System includes a Pre-Loaded Fathom-16 Steerable Guidewire into a Renegade HI-FLO Microcatheter. The Renegade HI-FLO Microcatheter incorporates a taper in its outside diameter along its length from the 3.0F (1.0 mm) proximal stiff region to the 2.8F (0.93 mm) flexible distal region. The I.D. of the microcatheter is 0.69 mm (0.027 in) minimally in the proximal and distal regions. The microcatheter lumen is able to accommodate steerable guidewires that are ≤ 0.47 mm (0.018 in) in diameter. The outer surface of the microcatheter's distal segment is coated with Hydro Pass™ Hydrophilic Coating. The Renegade HI-FLO Microcatheter has a radiopaque marker at the distal tip to facilitate fluoroscopic visualization. The distal tip of the microcatheter is steam shapeable. The proximal end of the microcatheter incorporates a standard luer with a rotating hemostatic valve (RHV) attached.

The FATHOM-16 Steerable Guidewire has a maximum diameter of 0.41 mm (0.016 in). The guidewire is compatible with existing microcatheters used in common endovascular procedures such as those for diagnosis and therapy in the peripheral vasculature. The distal portion of the device is coated with a hydrophilic polymer and the proximal portion is coated with polytetrafluoroethylene (PTFE). The distal 10 cm (3.9 in) of the guidewire tip is radiopaque. The distal 2 cm (0.79 in) is shapeable/reshapeable. The torque device is attached to the proximal end and functions as a steering mechanism. Rotation of this device facilitates guidewire placement into the appropriate vessel by precise directional manipulation of the guidewire tip.

COMPATIBILITY INFORMATION

Interface compatibility between any microcatheter, such as accessories and diagnostic, embolic or therapeutic agents for infusion, should be carefully considered before use. Consult tables below.

Table 1. Renegade HI-FLO Microcatheter Compatibilities

Guidewire	Guiding Catheter	EMBOLICS	
		Particles	Spherical
Max Dia. 0.47 mm (0.018 in)	Min. 0.96 mm (0.038 in) Guidewire compatible [(1.1 mm (0.042 in))]	≤ 710 Microns Results based on testing conducted using Contour™ PVA Embolic Particles	≤ 900 Microns Results based on testing conducted using Contour SE™ Microspheres

Table 2. FATHOM-16 Guidewire Compatibilities

Distal Outside Diameter	Proximal Outside Diameter	Compatible Microcatheters
0.41 mm (0.016 in) Maximum	0.41 mm (0.016 in) Maximum	Microcatheters with ID 0.53 mm – 0.69 mm (0.021-0.027 in) e.g. Renegade HI-FLO Microcatheter, Renegade™ Fiber Braided Microcatheter, Renegade™ STC-18 Microcatheter

Contents

This package contains one microcatheter, one guidewire, one torque device, one rotating hemostatic valve (RHV), one steam shaping mandrel, and one guidewire introducer.

INTENDED USE/INDICATIONS FOR USE

The Renegade HI-FLO FATHOM System is intended for peripheral vascular use. The FATHOM Guidewire can be used to selectively introduce and position the Renegade HI-FLO Microcatheter in the peripheral vasculature. The microcatheter can be used for controlled and selective infusion of diagnostic, embolic, or therapeutic materials into vessels.

CONTRAINDICATIONS

None Known.

ADVERSE EVENTS

The Adverse Events include, but are not limited to:

- Allergic reaction
- Death
- Embolism
- Hemorrhage/Hematoma
- Infection
- Pseudoaneurysm
- Stroke
- Vascular thrombosis

- Vessel occlusion
- Vessel spasm
- Vessel trauma (dissection, perforation, rupture)

WARNINGS

Never advance or withdraw an intravascular device against resistance until the cause of resistance is determined by fluoroscopy. Movement of the microcatheter or guidewire against resistance may result in damage or separation of the microcatheter or guidewire tip, or vessel perforation.

The infusion pressure with this microcatheter should not exceed 5520 kPa (800 psi). Infusion pressure in excess of this maximum may result in microcatheter rupture, possibly resulting in patient injury. If flow through the microcatheter becomes restricted, do not attempt to clear the microcatheter lumen by infusion. The static pressure with this microcatheter should not exceed 2070 kPa (300 psi). Static pressure in excess of this maximum may result in microcatheter rupture, possibly resulting in patient injury. Identify and resolve the cause of the blockage or replace the microcatheter with a new microcatheter before resuming infusion.

The Renegade™ HI-FLO™ FATHOM™ System is not intended for use in the coronary vasculature or the neurovasculature.

PRECAUTIONS

This device should be used only by physicians thoroughly trained in percutaneous, intravascular techniques and procedures.

Prior to use, carefully examine the unit to verify that the sterile package or product has not been damaged in the shipment.

Prior to a procedure, all equipment to be used for the procedure should be carefully examined to verify proper function and integrity.

Inspect the microcatheter and guidewire prior to use for any surface irregularities, bends or kinks. Any microcatheter or guidewire damage may decrease the desired performance characteristics.

Exercise care in handling of the microcatheter and guidewire during a procedure to reduce the possibility of accidental breakage, bending, or kinking. Do not use a microcatheter or guidewire that has been damaged.

When the microcatheter and/or guidewire are in the body, they should be manipulated only under fluoroscopy. Do not attempt to move the microcatheter and/or guidewire without observing the resultant tip response.

To avoid guidewire damage and possible shearing of plastic, do not withdraw or manipulate the guidewire through a metal needle cannula.

Extensive guidewire manipulation during lengthy procedures and the use of embolic agents may require the exchange of new microcatheters in place of the used microcatheter. Check labeled diameter of new microcatheter and verify compatibility with Table 2.

Because the microcatheter may be advanced into narrow subselective vasculature, repeatedly assure that the microcatheter has not been advanced so far as to interfere with its removal.

Excessive tightening of a hemostatic valve onto the microcatheter shaft or guidewire may result in damage to the microcatheter and/or abrasion of the coating of the wire.

Excessive tightening of the torque device onto the guidewire may result in abrasion of the guidewire PTFE coating.

Care should be taken when removing the Pre-Loaded System from the carrier tube as an exposed portion of the Pre-Loaded guidewire extends beyond the distal tip of the microcatheter.

It is recommended that a continuous saline flush be maintained between the guiding catheter and the microcatheter and between the microcatheter and the guidewire during the procedure. Flushing prevents contrast crystal formation and/or clotting on the guidewire and in the catheter lumen.

The recommended continuous flush set-up is shown in Figure 1.


Figure 1. Example of Continuous Flush Setup

The guidewire insertion tool, the RHV, the steam shaping mandrel, and the torque device are not intended to enter the body.

If other interventional devices are used with the Renegade HI-FLO FATHOM System, then refer to that product labeling for intended use, contraindications and potential complications associated with the use of that interventional device.

Do not use a product after the expiration date.

The Boston Scientific microcatheter, guidewire, torque device, rotating hemostatic valve (RHV), steam shaping mandrel, and a guidewire introducer are supplied STERILE and non-pyrogenic in unopened, undamaged packages. Verify that the sterility of the device has not been compromised by assuring package integrity has been maintained.

PREPARATIONS FOR USE

1. Connect a saline filled syringe to the carrier hoop luer.
2. Flush the catheter surface prior to removal from carrier hoop.
3. Remove angled clip and gently withdraw the microcatheter with the Pre-Loaded guidewire from the carrier tube. If resistance is felt during removal, re-flush with saline. Do not force removal from carrier tube. Once the catheter is hydrated, do not allow the product to dry. Do not reinsert product into the carrier tube.
4. Close RHV vent cap and remove proximal wire protector.
5. Flush microcatheter guidewire assembly with saline via RHV side arm.
6. Inspect guidewire and catheter prior to use to verify the microcatheter and guidewire are not damaged.

7. If desired, carefully shape the guidewire tip per the following instructions:
 - Make sure guidewire is hydrated.
 - Hold the guidewire introducer (included) or equivalent shaping tool in one hand and rest the guidewire tip across the wire introducer or shaping tool at a 90° angle.
 - Gently pinch the guidewire tip between the guidewire introducer or tool, and thumb.
 - Carefully pull the mandrel toward the guidewire distal tip to implement shape.
 - Repeat as needed to achieve desired shape.
 - Inspect the shaped tip carefully for any damage. If damaged, do not use the guidewire.
8. If desired, carefully shape the microcatheter tip per the following instructions:
 - Withdraw the guidewire so that the shaping mandrel can be inserted into the catheter's distal lumen.
 - Dip the distal end of the catheter into saline, bend the working length of the shaping mandrel to the desired catheter shape.
 - Shape catheter by holding mandrel/catheter assembly no closer than 25.4 mm (1 in) from a steam source for no more than 10 seconds. Multiple shaping is not recommended.
 - Cool tip in saline and remove the mandrel. Discard mandrel. Multiple shaping is not recommended.
9. Place the appropriate guiding catheter using standard technique. It is recommended that the Renegade™ HI-FLO™ Microcatheter be used with a guiding catheter that is 0.96 mm (0.038 in) guidewire compatible (with a minimum internal diameter of 1.1 mm (0.042 in)) and a sheath introducer.
10. Connect a rotating hemostatic valve (Tuohy-Borst type) to the guiding catheter luer per Figure 1 to continuously flush the guiding catheter with saline (a pressure pack is convenient for this purpose).

DIRECTIONS FOR USE

1. Align the distal tip of the guidewire to the distal tip of the microcatheter.
2. Carefully insert the microcatheter with the Pre-Loaded guidewire through the RHV into the guiding catheter per **Figure 1**. Tighten the RHV vent cap around the microcatheter to prevent backflow, but allowing some movement through the valve by the microcatheter.
3. If desired, the guidewire can be used with a torque device. To use the torque device, slip the device over the proximal end of the guidewire. When the torque device is in the desired location on the guidewire, tighten the cap to secure the torque device in place. The torque device may be repositioned or removed by loosening the cap and retightening the cap.
4. Set up a continuous flush using the RHV connected to the microcatheter (reference Figure 1). Ensure the RHV vent cap is closed and flush the microcatheter with Pre-Loaded guidewire.
5. Advance the guidewire and microcatheter to a selected vascular site by alternately advancing the guidewire and then tracking the microcatheter over the guidewire. This step is best accomplished by two people.
6. To rotate the tip of the guidewire to facilitate vessel selection, slowly and carefully rotate the torque device attached to the guidewire.
7. To infuse, completely remove the guidewire from the microcatheter. Connect a syringe with infusate to the microcatheter luer, and infuse as required.

8. When not in use during the procedure, wipe the guidewire with a heparinized saline and store in a saline bath. The guidewire should only be used for the same patient and procedure. For all agents, please refer to manufactures' instructions for use.

HOW SUPPLIED

Packaging is designed to maintain sterility according to the expiration date on the label. Do not use if package is opened or damaged. Do not use if labeling is incomplete or illegible.

Handling and Storage

Store in a cool, dry, dark place.

WARRANTY

Boston Scientific Corporation (BSC) warrants that reasonable care has been used in the design and manufacture of this instrument. **This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose.** Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond BSC's control directly affect the instrument and the results obtained from its use. BSC's obligation under this warranty is limited to the repair or replacement of this instrument and BSC shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this instrument. BSC neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. **BSC assumes no liability with respect to instruments reused, reprocessed or resterilized and makes no warranties, express or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instruments.**

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© 2015 Boston Scientific Corporation or its affiliates. All rights reserved.
Boston Scientific Limited
Boston Scientific (Australia) Pty Ltd