ANGIODYNAMICS®

VenaCure®
1470

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Version 2.0 April 2011 MAN/31/0075 EC

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- Foot Switch Assembly (supplied with laser):
AngioDynamics part number: AS1/A0/0002.
- Remote Interlock Lead (if required):
Lemo connector (supplied with laser): AngioDynamics part no. CON/51/0003
Ferrite sleeve (Farnell part no. 898-454).
Screened twisted pair cable 7/0.2 (Farnell part no. 140-457), maximum length 4m.

Note: The ferrite sleeve should be fitted to the remote interlock cable, at a maximum of 75 mm from the Lemo connector.

- Remote Interlock Bypass (supplied with laser):
AngioDynamics part number: AS1/A3/0024.
- IEC Mains Lead (supplied with laser):
AngioDynamics part numbers: CBL/02/0002, CBL/02/0040, CBL/02/0042, CBL/02/0046,
CBL/02/0051 or CBL/02/0063.

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mplications in extreme cases may include:

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Do not use coaxial gas/air coolant for non-contact fibres when there is a risk of air/gas embolism.

Do not use the laser close to large blood vessels or in highly vascularised areas, except when these are the target for the laser treatment.

LÂM SÀN**General Precautions**

Only operators who have been trained in the use of lasers and are thoroughly familiar with this Operator Manual should use the **VenaCure 1470**. The information provided in this section is not intended to be all-inclusive and it is not intended to replace operator training or experience. Please contact your AngioDynamics representative for training materials available on the use of this equipment.

Although it is difficult to specify the effect that the use of the diode laser will have in each therapeutic situation, it is possible to give a general overview as to what the clinician might expect when using the **VenaCure 1470**. The exact effect depends upon the chosen procedure and, especially when using the laser with manual settings, the Power setting, Pulse Duration, Pulse Interval, and the tissue type being treated.

The diode laser may cause tissue damage if improperly used. Precautions, such as careful assessment of the target tissue during treatment and the use of appropriate Power and Pulse Duration, should be taken. Use low Power and short Pulse Duration settings until fully familiar with the instrument's capabilities. Starting at low powers, the operator should note the effect on the tissue and increase Power, Pulse Duration or treatment time until the desired effect is obtained.


Specific parameters are not recommended, but are left to operator preference and best medical judgement dependent on the particular application.

The diode laser may not be effective for coagulation for severe haemorrhages. The operator must be prepared to control haemorrhages with strident, alternative non-laser techniques.

SAFETY LABELLING


Location of Safety Labelling

Safety labels for the **VenaCure 1470** are positioned as indicated below.

Product Identification Labelling

Product identification labelling is located on the rear of the **VenaCure 1470**.

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Guidance and manufacturer's declaration – electromagnetic immunity			
The VenaCure 1470 is intended for use in an electromagnetic environment specified below. The customer or the user of the VenaCure 1470 should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	+/- 6 kV contact +/- 8 kV air	+/- 6 kV contact +/- 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient / burst IEC 61000-4-4	+/- 2 kV for power supply lines +/- 1 kV for input/output lines	+/- 2 kV for power supply lines +/- 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	+/- 1 kV differential mode +/- 2 kV common mode	+/- 1 kV differential mode +/- 2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage Dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% U_T (>90% dip in U_T) for 0,5 cycle 40 % U_T (90% dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 sec	<5% U_T (>90% dip in U_T) for 0,5 cycle 40 % U_T (90% dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the VenaCure 1470 requires continued operation during mains interruptions, it is recommended that the VenaCure 1470 be powered from an uninterruptible power supply or battery.
Power frequency (50/60 Hz) IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE U_T is the a.c. mains voltage prior to application of the test level.			

Guidance and manufacturer's declaration – electromagnetic immunity

The **VenaCure 1470** Laser is intended for use in an electromagnetic environment specified below. The customer or the user of the **VenaCure 1470** Laser should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	<p>Portable and mobile RF communications equipment should be used no closer to any part of the VenaCure 1470, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance.</p> <p>$d = 1.2 \sqrt{P}$</p> <p>$d = 1.2 \sqrt{P}$ 80MHz to 800 MHz</p> <p>$d = 2.3 \sqrt{P}$ 800MHz to 2,5 GHz</p> <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).</p> <p>Field strengths from fixed RF transmitters as determined by an electronic site survey, ^a should be less than the compliance level in each frequency range. ^b</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,5 GHz	3 V/m	

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephone and land mobile radios, amateur radio, AM and FM radio broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the **VenaCure 1470** is used exceeds the applicable RF compliance level above, the **VenaCure 1470** should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the **VenaCure 1470**.

^b Over the frequency range 150 kHz to 800 MHz, field strengths should be less than 3 V/m.

Recommended separation distances between portable and mobile RF communications and the VenaCure 1470.

The **VenaCure 1470** is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the **VenaCure 1470** can help prevent electromagnetic interference by maintaining distance between portable and mobile RF communications equipment (transmitters) and the **VenaCure 1470** as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150kHz to 80 MHz $d = 1.2 \sqrt{P}$	80MHz to 800MHz $d = 1.2 \sqrt{P}$	800 MHz to 2,5 GHz $d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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It is essential that the exposed end of the optical fibre be kept clean to prevent damage to the VenaCure 1470 and optical fibre.

To remove the optical fibre connector, turn the gripper anti-clockwise until fully unscrewed and disconnect from the laser aperture. Dispose of the optical fibre according to institution policy. Alternatively, if the optical fibre is permitted for multiple uses, immediately fit a protective cap over the end of the optical fibre to protect the optical surface from contamination.

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¹ The system will read the fibre's data if the power is on and a FRS fibre is connected or if the **VenaCure 1470/30** is powered up with a FRS fibre already connected.

² The number of available "uses" is only decreased once per procedure and only after the **VenaCure 1470/30** has been fired. Subsequent firing of the **VenaCure 1470/30** during a continuous procedure will not result in the number of available "uses" being further decreased.

³ Attaching an incompatible FRS fibre will result in no procedures being available for selection in the main menu.

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The method of adjusting settings will vary depending on the procedure selected. Please refer to section 4 entitled *AngioDynamics Procedures* for full details of how to use specific procedures.

⁴ The list of procedures shown on this menu will only include those that are installed in the laser. AngioDynamics will, over the lifetime of this system, release new procedures that may be compatible with currently available fibres. Please check with your AngioDynamics representative for the latest information regarding available procedures. Further information on the currently available procedures can be found in section 4.

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Description	AngioDynamics part number
EVL TM Test fibre	AS/0604
Manual Test fibre	AS/0695
Emergency Override Reset Device	AS/0598
Laser safety glasses (premium quality)	SE/0005
Laser safety glasses (standard quality)	SE/0006
Laser warning sign	PP/1580
Memory Card	S10/06/0086
USB Memory Card Reader/Writer	S10/06/0087

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
Theo bảng sau:


	EVLT™ Continuous	EVLT™ Manual
Mode(s)	Continuous	Continuous or Repeat Pulse
Power (W) 	6W (4 W to 10 W)	1 W to 12 W
Pulse Duration (s) 	Not available	0.1 s to 2 s
Interval (s) 	Not available	0.1 s to 2 s

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




	PVAK Continuous	PVAK Manual
Mode(s)	Continuous	Continuous, Pulsed or Repeat Pulse
Power (W) 	6 W (4W to 10W)	1 W to 12 W
Pulse Duration (s) 	Not available	0.1 s to 5 s
Interval (s) 	Not available	0.1 s to 5 s

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SPECIFICATIONS

VenaCure 1470	
Laser Type	Diode laser, CW
Centre wavelength	1470 nm
Spectral bandwidth	90% of optical power contained within 1470 ± 20 nm
Delivery fibre	Minimum 400 µm diameter, 0.37 Numerical Aperture (NA)
Fibre connector	Optical SMA-905 to MIL STD
Maximum Power	12 W to the laser output port
Minimum Power	1.0 W to the laser output port
Set power accuracy	Better than ± 10% of displayed power
Output power stability	± 5% maximum
Aiming beam	Red Class 3R diode laser, ≤ 5 mW at port, wavelength 635-660 nm
Pulse duration	100 ms – 9900 ms
	Continuous Wave (CW)
Pulse Interval	100 ms – 1000 ms
Fibres	VenaCure EVLT™ kits only
Fibre Recognition System (FRS)	
Operating Frequency	125 kHz ± 1 kHz
Operating Range	Maximum of 40 mm (1.6") from the laser output port
RFID Modulation Type	A1D
Standards	FCC 47 Part 15c, Canadian RSS-210 EN 300 330-1, EN 300 330-2, EN 301 489-1, EN 301 489-3
Cooling	By ambient air, with fan assistance
Power Supply	100 – 240 V AC, 50 – 60 Hz, 500 VA max.
Dimensions (H x W x D)	150 x 370 x 325 mm (± 5 mm) (5.9" x 14.6" x 12.8" (± 0.2"))
Weight	12 kg max. (26.4 lb)
Operating Temperature	10°C to 30°C (50°F to 104°F)
Operating Humidity	Up to 75% relative humidity, non-condensing
Operating Pressure	Normal atmospheric pressure
Storage Temperature	0°C to 55°C (32°F to 130°F)
Storage Humidity	Up to 90% relative humidity, including condensing
Storage Pressure	500 hPa to 1014 hPa
Safety Standards	UL 60601-1, 21 CFR 1040.10, 1040.11, EN 60601-1, 60601-1-2, 60601-2-22, 60825-1,

CLASSIFICATION

Type of protection from electric shock	Class I
Degree of protection from electric shock	Type B
Degree of protection against ingress of water	IP21

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FRS Troubleshooting

The AngioDynamics FRS system has been extensively tested. In the unlikely event of a problem being experienced please perform the following checks before contacting your local AngioDynamics representative for further advice.

- “Invalid Fibre” is displayed even with a new unused FRS system fibre**
 - (1) Disconnect and then reconnect the fibre to the **VenaCure 1470**. This will cause the unit to try and read the fibre’s data again.
 - (2) Switch the **VenaCure 1470** off and on at the power inlet.
 - (3) Ensure that the **VenaCure 1470** is at least 2 meters away from any other electrical or electronic equipment that might interfere with the FRS system, such as computers or other electronic or medical equipment.

I need to use a non-FRS system fibre The **VenaCure 1470** will normally only operate with AngioDynamics FRS fibres. However, if the ‘u’ symbol shown right is displayed in the bottom left corner of the screen, then other non-FRS fibres may also be used. Please contact your AngioDynamics representative for further details of this feature.

The Emergency Override option is not available After it has been used once, the Emergency Override option must be reset before it can be used again, as described in ‘Resetting the Emergency Override’ above. Contact your AngioDynamics representative to obtain a spare reset device.

The fibre is still within its sterility date but the laser shows it as expired Check that the date of the internal clock in the **VenaCure 1470** is set correctly. The procedure for setting the clock is described in Section 3 - Operator Instructions.

ACCESSORIES***Optical Fibres***

The **VenaCure 1470** has an output connector for optical fibres with standard SMA-905 connector. Only AngioDynamics labelled fibres should be used. A list of fibres available for use with the **VenaCure 1470** can be obtained from your AngioDynamics representative.

Test Fibre

The test fibre is a special type of FRS fibre used during testing and evaluation only.

Emergency Override Reset Device

The Emergency Override Reset Device consists of a special FRS gripper, supplied without a fibre. Instructions on when and how to use this device are in section 3.

AngioDynamics does not recommend the use of third party fibres as their quality and efficacy cannot be guaranteed. Any damage caused to your **VenaCure 1470** by using a fibre not supplied by AngioDynamics may not be covered under the AngioDynamics warranty.

STERILIZATION OF OPTICAL FIBRES

Optical fibres are provided sterile as a disposable, single-use product.

DO NOT RE-STERILIZE THE FIBRES.

DO NOT RE-USE THE FIBRES.

USE ONLY ANGIODYNAMICS LABELLED OR ANGIODYNAMICS APPROVED FIBRES.

Failure to observe this could invalidate the Laser Warranty.

Intra-operative cleaning

If the tip accumulates debris, turn the laser to the STANDBY mode and then carefully wipe the tip clean with a wet sponge/swab.

Fibre Disposal

After use, the single-use optical fibres should be disposed of in accordance with local regulations regarding disposal of contaminated waste.

AngioDynamics labelled optical fibres have undergone stringent evaluation and testing to ensure that they are of the highest quality and that they operate safely, effectively and efficiently with AngioDynamics lasers.

The exact alignment of the interface between the laser aperture and the SMA-905 connector is critical. Misalignment (as may occur with non-approved fibres) can result in damage to the laser and poor delivery of laser energy to the patient.

FUSE REPLACEMENT

Spare power fuses are supplied with the **VenaCure 1470**. Further spares can be obtained from AngioDynamics. They can be replaced as follows:

1. Disconnect the **VenaCure 1470** from the AC supply.
2. Use a small flat-bladed screwdriver to release the fuse compartment from the mains inlet on the rear panel.
3. Remove the two fuses from the holder and replace with new ones of the same type and rating: T6.3A H 250V. Fuses with a different rating or specification must not be used.

DISPOSAL

At the end its life the **VenaCure 1470** should be disposed of according to national environmental requirements or be returned to AngioDynamics.

SOFTWARE UPDATES

From time to time AngioDynamics may issue new procedures, upgrades and feature enhancements for the **VenaCure 1470**. As a registered owner you will be notified of these when they become available.

SERVICING

The **VenaCure 1470** does not require regular servicing or maintenance, with the exception of annual checks of the calibration of the output power of the laser. If either of these checks is due, a 'spanner' symbol will be displayed on the screen at start-up. This symbol will be removed after the checks have been completed.

LASER POWER OUTPUT

Measuring the Laser Power Output

The LASER SAFETY OFFICER or suitably trained service personnel should check the output power of the **VenaCure 1470** at least annually from the date of installation, by following the procedure described below.

- Equipment Required**
- A sampling power meter or an independent energy (integrated power) meter of known calibration
 - A bare ended optical fibre
 - Laser unit to be tested

Procedure

Calibrate the fibre.
 Connect the fibre to the Laser unit output port and present the distal end of the fibre to the external power meter.

1. Record the Laser unit's actual and displayed outputs at various different power/energy settings e.g. 5W, 10W etc.
 Calculate the percentage difference between the displayed and the actual power/energy output as taken from the external power meter.

Adjusting Laser Power Output

2. If calculated disparity exceeds $\pm 20\%$, contact AngioDynamics.
 Power output adjustments can only be made by suitably trained AngioDynamics service personnel.
 For regulatory purposes, the method for carrying out these adjustments is described below.
 Please contact your AngioDynamics representative for further advice.

Adjusting the Laser Power Output

THE FOLLOWING INFORMATION IS PROVIDED FOR REGULATORY PURPOSES

The procedure below will explain how to adjust the power output of the laser.

- Equipment Required**
- A sampling power meter or an independent energy (integrated power) meter of known calibration
 - A calibrated bare ended optical fibre
 - Laser unit to be tested

Procedure

Connect the laser to a PC via a null-modem serial cable between the port on the rear panel and a serial port on the PC.
 Connect the calibrated fibre from the laser unit output port to the external power meter.
 Turn on the laser.
 Load and run the AngioDynamics Engineer Interface program on the PC. From the **Connection** menu, select **Connect to Laser**, choose the appropriate Comms Port and press **Connect**.
 From the **Connection** menu select **Enter EI** and enter the password.
 Select the **Non-volatile** data tab.
 Diode calibration is stored in both the Monitor MCU and in the laser module itself. The values in these locations are shown on the screen – if there are any discrepancies then an error message is displayed.
 The laser is calibrated at three different power settings: 2W, 5W and 12W, by entering the required diode feedback value (in mV) into the appropriate boxes on the form. Always change the values for both the **Monitor params** and **Diode data** at the same time, then press **Update** to reconfigure the laser.

Start with the 2W setting by entering a default value (no more than 150) into the two appropriate boxes and pressing **Update**.

Calculate the expected power at the distal end of the fibre, i.e. measured at the power meter, from the required port power and the known efficiency of the fibre.

Fire the laser and record the power measured by the power meter.

Adjust the value of the feedback according to this result and enter new values until the correct power is observed when the laser is fired.

Repeat this process for the 5W and 12W settings, starting with default values of 400 and 1000 respectively.

Finally, check that the laser outputs the correct power at all three settings.

TECHNICAL DESCRIPTION

The **VenaCure 1470** laser contains no user-serviceable components. In the event that repair or service is required please contact your local AngioDynamics representative. More detailed service instructions, including schematic diagrams, are available only to suitably-qualified and trained technical personnel.

The **VenaCure 1470** laser is built in a modular fashion, enabling ease of test, assembly and service. The system has been designed so that any one module can be replaced with no performance effect on the other modules or the product as a whole.

The modules present in a complete unit are:

- Laser Module
- Power Supply / Laser & TEC Driver
- Monitor & Control PCA
- User Interface PCA

Laser Power Control Overview

Optical power control is achieved using a monitor photodiode measuring the output of the laser diode. This is used to derive a control signal proportional to the total power output of the laser unit. The feedback signal in this control loop is monitored for errors from the expected value.

In addition to this, an over-current trip circuit will operate to disable the laser driver rapidly should an overcurrent / overpower situation be detected. Diode current is monitored while firing to ensure that it is within acceptable limits.

Laser Module

The Laser Module consists of a laser diode, TECs, optics and an electronic control module. These components are all mounted on a stable metal platform which doubles as a heatsink. Two fans are mounted at one end of the heatsink to provide air-assisted cooling. The laser and TECs are located inside a hermetically-sealed enclosure.

Semi-transmissive mirrors in the optics path allow for the addition of the visible laser diode and for output power to be measured by means of a photodiode. The laser output port (SMA) incorporates a thermistor, to measure the temperature of the port close to the fibre connection and two micro-switches, which detect when a fibre has been correctly connected to the port.

PSU / Laser & TEC Driver

The power supply, laser driver and TEC driver are contained in their own enclosure. The power supply accepts an input from 100V to 240V and provides auxiliary power rails for the system electronics. The drivers are controlled via an interface with the Monitor & Control PCA.

Monitor & Control PCA

The Monitor & Control PCA consists of the laser power control loop, Control and Monitor microcontrollers, control logic, external interlock interfaces, power supply monitoring, cooling fan control and the RS232 serial port connection. This PCA connects to the Laser Module, the User Interface and the Power Supply.

Two separate microcontrollers are used to independently Control and Monitor the laser diode. Instructions and information are passed from these two devices to the User Microcontroller via a dedicated RS232 serial interface.

Monitoring of the laser diode and peripheral functions such as temperature is performed by the Monitor MCU. The values obtained are also sent to the User Microprocessor. The Monitor MCU additionally

monitors the state of many of the Control MCU outputs to check that their status is correct for the operating mode.

The Footswitch and Remote Interlock interface connectors are located on the rear panel. All signals on these connectors are filtered for EMC and protected against ESD. The Remote Interlock connector also includes an electrically-isolated 'READY OUT' output.

The two cooling fans attached to the heatsink may be run at four different speeds, depending on the temperature of the heatsink, the power of the laser output and the magnitude of the TEC drive level.

The signals for the RS232 serial port come direct from the User Microcontroller. The pins of the 9-pin D-type connector are electrically isolated from the rest of the circuit using opto-couplers and a transformer.

Ambient temperature is monitored using a digital temperature sensor, which is read by the User Microcontroller over the I²C bus.

User Interface PCA

The User Interface PCA consists of the User Microcontroller, display and all the controls and indicators required by an operator to use the product.

The User Microcontroller is from the Renesas H8 family. It is supplemented by 256Kb of SRAM memory, 4Mb of Flash program memory and 256Kb of serial FRAM memory for non-volatile data storage. A Real Time Clock (RTC) function provides the time and date features for the user interface. This microcontroller system interfaces to the rest of the product via a Serial Peripheral Interface (SPI) bus, Inter-Integrated Circuit (I²C) bus and RS232 serial interface.

The user interface display is a Vacuum Fluorescent Display (VFD) of 128 x 64 pixels. This is controlled directly by the User Microcontroller. The display is filtered to mid-blue.

The control of the level of TEC drive is a software function within the User Microcontroller. It operates in a feedback loop with the objective of keeping the laser diode at a fixed temperature.

A memory card interface allows a MMC or SD memory card to be inserted into an aperture on the front panel. This card is used to transfer usage data from the laser to the user's PC. It can also be used to transfer updated software and new procedure information into the laser.

Glossary

ADC	Analogue to Digital Converter
DAC	Digital to Analogue Converter
EMC	Electro-Magnetic Compatibility
EEPROM	Electrically Erasable Programmable Read Only Memory
ESD	Electro-Static Discharge
FRS	Fibre Recognition System
I ² C	Inter-Integrated Circuit
I/O	Input/Output
Kb	kilobits
LED	Light Emitting Diode
MCU	Microprocessor Control Unit
ms	milliseconds
MMC	MultiMediaCard
OPCM	Optical Power Calibration Meter
PCA	Printed Circuit Assembly
PCB	Printed Circuit Board
PSU	Power Supply Unit
PWM	Pulse Width Modulation
RAM	Random Access Memory
RTC	Real Time Clock
SD	Secure Digital memory card
SFC	Single Fault Condition
SPI	Serial Peripheral Interface
TTL	Transistor-Transistor Logic
TEC	Thermo Electric Cooling
UART	Universal Asynchronous Receiver/Transmitter
V	Volts
VFD	Vacuum Fluorescent Display

SECTION 6 GUARANTEE

MANUFACTURER'S GUARANTEE POLICY

AngioDynamics guarantees the **VenaCure 1470** against defects in materials and workmanship for a period of 1 year. The guarantee period begins on the date of purchase.

To enable timely registration of the guarantee, the owner/purchaser must complete and return the Guarantee registration form located in the guarantee booklet within 28 days of purchase.

The following items are expressly excluded from this Guarantee:

- Safety Eyewear
- All optical fibres and accessories
- Maintenance Instruments
- All other accessories supplied by AngioDynamics

Any attempt to repair, adjust or modify the system beyond those procedures described in the Operator Manual by any person not authorized by AngioDynamics, will invalidate the Guarantee.

**GUARANTEE
CLAIMS**

To make a guarantee claim the purchaser shall, promptly following discovery of the basis of claim, contact their AngioDynamics distributor in the first instance or AngioDynamics UK Ltd. in writing, by telephone, fax or Email at the following address:

ANGIODYNAMICS®

AngioDynamics UK Limited
Building 2000, Beach Drive
IQ Cambridge, Waterbeach
Cambridge CB25 9TE
United Kingdom

Tel: +44 1223 729372

Tel: 0800 013 0859 (UK only - free phone)

Fax: +44 1223 729320

Email: customerservice@angiodynamics.com

AngioDynamics is a Registered Trademark of AngioDynamics Inc.

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© 2011 AngioDynamics UK Ltd.
VenaCure, VenaCure EVLT, NeverTouch and EVLT are Registered Trademarks of AngioDynamics Inc.