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angiodynamics

INDICATIONS:

The AngioDynamics*, Inc. VenaCure EVLT* NeverTouch Direct* Procedure Kits are indicated for endovascular coagulation of the Great Saphenous Vein (GSV) in patients with superficial vein reflux, for the treatment of varicose veins and varicosities associated with superficial reflux of the Great Saphenous Vein (GSV), and for the treatment of incompetence and reflux of superficial veins of the lower extremity.

This product should be used only with lasers cleared for use in the treatment of varicose veins, varicosities with superficial reflux of the GSV, and in the treatment of incompetent refluxing veins in the superficial venous system in the lower limbs.

EXCLUSIONS:

- Patients with thrombus in the vein segment to be treated
- Patients with an aneurysmal section in the vein segment to be treated
- Patients with peripheral artery disease as determined by an Ankle-Brachial Index <0.9

POTENTIAL COMPLICATIONS:

The potential complications include but are not limited to the following: vessel perforation, thrombosis, pulmonary embolism, phlebitis, hematoma, infection, skin pigmentation alteration, neovascularization, paresthesia due to thermal damage of adjacent sensory nerves, anesthetic tumescence, non-target irradiation, hemorrhage, necrosis, skin burns and pain.

PRECAUTION:

- Prior to and during use, avoid damaging the fiber by striking, stressing, or excessive bending. Do not coil the fiber tighter than a diameter of 16cm. Clinical safety and effectiveness data is not available for other fiber tip designs and diameters.
- Prior to and during use, avoid bending the introducer sheath and dilator as this can cause kinks and damage.
- All persons in the treatment room must be wearing protective glasses with the proper rating for the wavelength being used.

WARNINGS:

- Treatment of a vein located close to the skin surface may result in skin burn.
- Paresthesia may occur from thermal damage to adjacent sensory nerves.
- Tissue not targeted for treatment must be protected from injury by direct and reflected laser energy. Use appropriate eye and protective wear for both patient and operating personnel.
- Reuse of single-use devices creates a potential risk of patient or user infections. Contamination of the device may lead to injury, illness, or death of the patient.
- Reprocessing may compromise the integrity of the device and/or lead to device failure.

INSTRUMENTS REQUIRED:

- Laser
- Safety eyewear

CONTENTS:

- NeverTouch Direct Fiber with SMA 905 connector
- 4 French, 10cm Introducer
- 0.018 in, 45cm Guidewire
- 21 Gauge, 7cm Needle

FIBER REMOVAL FROM DIE CARD:

NOTE: Orient the die card such that the AngioDynamics logo is legible.

1. Use aseptic technique to remove all components from pouches.

A. To remove fiber tip end first:

- 2a. Advance the fiber tip into the ship tube slightly to facilitate removal from die card tab.
- 3a. Unclip the fiber from the bottom right inside tab of the card.
- 4a. Pull the fiber tip out of the ship tube.
- 5a. Uncoil the fiber in a counterclockwise manner, pulling the fiber vertically from the card.
- 6a. Unclip the SMA from the raised tabs.

B. To remove fiber SMA end first:

- 2b. Unclip the SMA from the raised tabs of the card.
- 3b. Uncoil the fiber in a clockwise manner, pulling vertically from the card.
- 4b. Pull the fiber tip from the ship tube.

PROCEDURE:

Carefully read all instructions and observe all Warnings prior to performing the procedure. Patient complications may occur if this is not done. This procedure can be performed in the physicians' office or as an outpatient treatment under local anesthesia and should be performed by a qualified physician who has received training in these techniques.

1. Carry out a physical examination.
2. Review ultrasound vein diameter and depth measurements.
3. Mark the skin with a surgical pen to show the length of vein to be treated.
4. Set the laser up according to the provided Diode Laser Operator Manual.
5. Place the patient in reverse Trendelenburg position to engorge and distend the targeted treatment vein.
6. Using standard Seldinger technique, gain access to the desired vein with the needle and corresponding guidewire.
7. Advance the introducer/dilator over the wire.
8. Remove the dilator and guidewire.
9. Properly attach the sterile, single use disposable NeverTouch Direct laser fiber by removing the protective cap from the fiber fitting and completely screwing the fiber fitting clockwise to the laser's SMA 905 connector until tight and secure.
10. Advance the NeverTouch Direct laser fiber through the introducer sheath to the treatment location. If treating the Great Saphenous Vein, position the fiber tip so it is 1-2cm below the Sapheno-Femoral Junction. Verify fiber tip position using ultrasound guidance.
11. Administer local tumescent anesthesia under ultrasound guidance to ensure delivery of just enough anesthesia to achieve thermal protection.
12. Set the laser to continuous mode and adjust to desired power.
(810nm, 980nm Lasers: 10-14 Watts, 1470nm Lasers: 5-7 Watts)
13. Put on laser safety glasses appropriate for the applicable wavelength.
14. Place the laser in the "Enable" mode.
15. Slide the introducer sheath out of the vessel and back along the fiber.
NOTE: If the sheath is left inside the vessel during the procedure, upon visualization of the 16cm warning mark, begin removing the sheath from the vessel in conjunction with the fiber.
NOTE: If use with the 65 cm Trè-Sheath* introducer is required in place of the introducer, align the back end of the Trè-Sheath hub with the center of the 72cm locating mark.
NOTE: If a Trè-Sheath is used, confirm the fiber tip protrudes at least 2cm past the tip of the Trè-Sheath prior to lasing.
16. Verify fiber tip position using ultrasound guidance.
17. Activate the laser by depressing the foot pedal while withdrawing the fiber (and Trè-Sheath introducer if applicable), at a rate that adequately delivers desired laser energy per centimeter.
(810nm, 980nm Lasers: 50-80 Joules per cm)
(1470nm Lasers: 30-50 Joules per cm)
Do not apply direct external hand pressure or force over the fiber tip during energy activation.
18. Cease operating the laser by removing foot from the foot pedal when the fiber tip is 2 – 3 cm from the access site as indicated by markers on the shaft of the fiber (or Trè-Sheath introducer if applicable). The NeverTouch Direct fiber tip is 4 cm from the access site when the white exit marks begin and 3cm from the access site when the white exit marks terminate.

The procedure is considered complete when the desired length of the targeted treatment vein is treated.

POST PROCEDURE:

Perform a second ultrasound to confirm that the vein is closed and that there is no more blood flow.

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AngioDynamics, VenaCure EVLT, Trè-Sheath and NeverTouch Direct are trademarks and/or registered trademarks of AngioDynamics Inc., an affiliate or a subsidiary. Covered by or for use under one or more of the following U.S. Patents: 7,559,329; 7,483,457; 7,907,643; 8,425,501; 8,864,754; 8,864,755; other U.S. patents and foreign counterparts pending.
©AngioDynamics, Inc. 2019 all rights reserved.