# **Transend™**

Guidewire with  
ICE™ Hydrophilic Coating

<b>Directions for Use</b>	<b>2</b>
<b>Instrucciones de uso</b>	<b>5</b>
<b>Mode d'emploi</b>	<b>9</b>
<b>Gebrauchsanweisung</b>	<b>12</b>
<b>Istruzioni per l'uso</b>	<b>16</b>
<b>Gebruiksaanwijzing</b>	<b>19</b>
<b>Instruções de Utilização</b>	<b>23</b>

90959072-01

2014-08

# Transend™

## Guidewire with ICE™ Hydrophilic Coating

### **Rx ONLY**

**Caution:** Federal Law (USA) restricts this device to sale by or on the order of a physician.

---

#### **WARNING**

Contents supplied STERILE using an ethylene oxide (EO) process. Do not use if sterile barrier is damaged. If damage is found, call your Boston Scientific representative.

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

---

#### **DEVICE DESCRIPTION**

The Transend Guidewire is a steerable guidewire available in a nominal diameter of 0.014 in (0.37 mm) or 0.018 in (0.47 mm). The Transend 14 Guidewire is available in nominal lengths of 135, 165 and 190 cm and Transend 18 Guidewire is available in nominal lengths of 135 and 165 cm. The distal 38 cm of the guidewire contains a core wire that is jacketed with a radiopaque polymer and is coated with ICE hydrophilic coating. The distal 2 cm of the guidewire is shapeable. The proximal portion of the guidewire is coated with Polytetrafluoroethylene (PTFE). A torque device (pin vise) is included with the wire to facilitate directional manipulation of the guidewire. A guidewire introducer is included with the wire to facilitate the introduction of the guidewire into catheter hubs or hemostatic valves.

#### **Contents**

- One (1) Insertion Tool
- One (1) Torque Device
- One (1) Guidewire

#### **INTENDED USE/INDICATIONS FOR USE**

The Transend Guidewire is intended for general intravascular use, including the peripheral vasculature. The wire can be torqued to facilitate the selective placement of diagnostic or therapeutic catheters.

#### **CONTRAINDICATIONS**

This device is not intended for use in coronary arteries.

#### **PRECAUTIONS**

- This device should be used only by physicians thoroughly trained in percutaneous, intravascular techniques and procedures.
- Prior to a procedure, all equipment to be used for the procedure should be carefully examined to verify proper function and integrity.
- Inspect guidewire prior to use for any surface irregularities and bends or kinks. Damaged and/or irregular guidewires should not be used.

- Exercise care in handling of the guidewire during a procedure to reduce the possibility of accidental breakage, bending, or kinking. Do not use a guidewire that has been damaged.
- To avoid guidewire damage and possible shearing of plastic, do not withdraw or manipulate the wire through a metal needle cannula.
- When the guidewire is in the body, it should be manipulated only under fluoroscopy.
- Do not attempt to move the wire without observing the resultant tip response.
- Never advance or withdraw an intravascular device against resistance until the cause of the resistance is determined by fluoroscopy. Excessive force against resistance may result in separation of the guidewire tip, damage to the catheter, or vessel perforation.
- Due to the variations of certain catheter tip inner diameters, abrasion of the hydrophilic coating may occur during manipulation. If any resistance is felt during introduction of the catheter, use of a different catheter may be warranted.
- Neither the guidewire insertion tool, nor the torque device are intended to enter the body.
- Excessive tightening of the torque device onto the wire may result in abrasion of the coating on the wire.
- The Boston Scientific guidewire, the guidewire insertion tool, and the torque device are supplied STERILE and non-pyrogenic in unopened, undamaged packages. Verify that the sterility of the device has not been compromised by assuring package integrity has been maintained.
- If the integrity of the packaging has been compromised, **do not use or attempt to resterilize**. Contact your local Boston Scientific representative
- Check labeled diameter of diagnostic or therapeutic catheter and verify compatibility with the guidewire outer diameter prior to use.

#### **ADVERSE EVENTS**

Complications attributed to guidewire applications are the following:

- Procedural related complications including but not limited to:
  - Vessel trauma
  - Vessel damage
  - Air embolism, thromboembolism
  - Post embolization syndrome (abdominal pain, fever, and nausea/vomiting)
  - Hematoma at the puncture site
  - Infection
  - Perforation of the vessel
  - Vessel spasm
  - Hemorrhage
  - Vascular thrombosis
  - Death
  - Bleeding
- Failed treatment
- Inability to position guidewire
- Damage to catheter
- Excessive force against resistance may result in separation of the guidewire tip

## HOW SUPPLIED

### Handling and Storage

Store in a cool, dry, dark place. Do not use if package is opened and/or damaged. Use the device prior to the "Use By" date noted on the product label. Do not use if labeling is incomplete or illegible.

### PREPARATION FOR USE

Before removing the guidewire from the carrier tube, flush the carrier tube with heparinized saline to hydrate the hydrophilic segment of the wire. Repeat injection if difficult removal of the wire occurs. Remove the guidewire carefully from the carrier tube to reduce the possibility of damage to the distal tip.

### DIRECTIONS FOR USE

1. Prior to inserting the guidewire into a catheter, flush the catheter with heparinized saline. This will prime the catheter and also provide smooth movement of the guidewire within the catheter.
2. Gently remove the guidewire from the carrier tube and inspect the wire prior to use to verify that it is undamaged.
3. If desired, gently shape the guidewire flexible tip according to standard practices.
4. Carefully insert and advance the guidewire into the catheter. The guidewire insertion tool provided may be used to facilitate the introduction of guidewires into catheter hubs or hemostatic valves.
5. If desired, the guidewire can be used with the torque device provided. To use the torque device, slip the device over the proximal end of the guidewire. When the torque device is in the desired location on the guidewire, tighten the cap to secure the torque device in place. The torque device may be repositioned or removed by loosening the cap and retightening the cap.
6. To rotate the tip of the guidewire to facilitate vessel selection, slowly and carefully rotate the torque device attached to the guidewire.
7. When not in use during the procedure, wipe the guidewire with heparinized saline and store in the carrier tube in a saline bath. This guidewire should only be used for the same patient and procedure.

### WARRANTY

Boston Scientific Corporation (BSC) warrants that reasonable care has been used in the design and manufacture of this instrument. **This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose.** Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond BSC's control directly affect the instrument and the results obtained from its use. BSC's obligation under this warranty is limited to the repair or replacement of this instrument and BSC shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this instrument. BSC neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. **BSC assumes no liability with respect to instruments reused, reprocessed or resterilized and makes no warranties, express or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instruments.**

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All rights reserved.
Boston Scientific Limited
Boston Scientific (Australia) Pty Ltd
© 2014 Boston Scientific Corporation or its affiliates.