
CASCADETM

Non-Occlusive
Remodeling Net

INSTRUCTIONS FOR USE

CASCADE™

Non-Occlusive Remodeling Net

READ THIS PAMPHLET CAREFULLY IN ITS ENTIRETY BEFORE USING THE DEVICE!

This document is intended to guide the operator in the use of this product. It is not intended to serve as a reference for endovascular coiling or other endovascular techniques. Refer to all warnings, cautions, contraindications, and instructions prior to use. Please save this pamphlet for further reference.

DESCRIPTION

The Cascade Non-Occlusive Remodeling Net is a sterile, non-pyrogenic and flexible device that provides mechanical support during coil embolization of intracranial aneurysm, by preventing coils from herniating into the parent vessel, while allowing continuous blood flow in the parent vessel. The device is designed with a braided net that can be adjusted to the vessel diameter, via the actuation handle (Figure 1). The actuation handle has a free mode feature, allowing for continuous adjustments, and an auto-lock mode, which allows for adjustments in set increments. The device is inserted and deployed under fluoroscopic visualization and is compatible with catheters and micro-catheters of 0.017" inner diameter (Cascade 17), 0.021" inner diameter (Cascade M, L, and Agile models), or larger. Once positioned, the Cascade can be expanded in full compliance with the parent vessel's shape and diameter, assisting in positioning the coiling microcatheter and coil deployment as needed. The Cascade is available in several sizes and models, providing coverage for vessels ranging between 1.5 – 6 mm in diameter, and aneurysms with neck width of up to 10mm.

INDICATIONS FOR USE

The Cascade is intended to provide temporary assistance for coil embolization of intracranial aneurysms.

CONTRAINDICATIONS

- > Patients with a known hyper-sensitivity to nickel-titanium.

CAUTIONS

- > The Cascade should only be used by physicians trained in interventional neuroradiology and treatment of intracranial aneurysms.

- > Use the Cascade in conjunction with fluoroscopic visualization agents throughout the entire procedure, in accordance with standard medical practice.
- > Use the Cascade in conjunction with anticoagulant and anti-platelet agents, in accordance with standard medical practice.
- > Check all connections to make sure that no air enters the guide catheter or microcatheter during continuous flushing.
- > Keep the Cascade device hydrated throughout the procedure.

WARNINGS

- > Before use, refer to the complete warnings, contraindications and instructions that follow.
- > When using the instrument, sound medical practice must always be followed by the user.
- > The Cascade is a delicate instrument and must be handled carefully. Before use and when possible, during the procedure, carefully inspect device integrity. Do not use a device that shows signs of damage. Damage may prevent the device from functioning and may cause complications.
- > Use a new microcatheter for each new Cascade device.
- > Do not rotate or torque the device to avoid injury or device damage.
- > Do not attempt to re-shape the Cascade's tip, as it might damage the device.
- > The Cascade is provided sterile and should only be used in a sterile environment
- > Visually inspect the sterile package for breaches of packaging integrity prior to use. Verify that it has not been damaged during shipping. Do not use if package is open or damaged.
- > Do not re-sterilize and/or reuse. Reprocessing and reuse increase the risk of infection and compromise device integrity and performance.
- > Do not autoclave or expose to extreme temperatures as it may damage the device.
- > Do not expose the Cascade to any detergents.
- > Do not use beyond "Use By" date.

COMPLICATIONS

Possible risks include, but are not limited to: air embolism, hematoma or hemorrhage; infection; distal embolization; vessel spasm, thrombosis, dissection, or perforation; emboli; acute occlusion; ischemia; neurological deficits including stroke; death.

PREPARATION PROCEDURE

1. Remove the Cascade from the package hoop.
2. Carefully unsheathe the braided net from the loading tube by retrieving the loading tube backwards.
3. Check that the Cascade functions properly by actuating the handle and verifying corresponding diameter changes in the braided net (see figure 2).
4. Verify that both modes of handle operation are effective (Free Mode and Auto Lock). Switch between modes by rotating the top part of the handle by 90° (see Figure 3).
5. Gradually decrease the diameter of the braided net while carefully re-sheathing with the loading tube.


Figure 1 - Cascade parts

Table 1 - Device properties and sizes

Model	Effective Braid Diameter (A)	Maximal Braid Length (B)	Usable Length (C)	Tip to Flurosaver Marking Distance (D)	Recommended Vessel Diameter	Min. Micro-catheter ID
Cascade 17	0.4 – 3 (mm)	18 (mm)	191 (cm)	150 (cm)	1.5 – 3 (mm)	0.017 (inch)
Cascade M	0.5 – 4 (mm)	37 (mm)	193 (cm)	150 (cm)	2 – 4 (mm)	0.021 (inch)
Cascade M Agile	0.5 – 4 (mm)	21 (mm)	191 (cm)	150 (cm)	2 – 4 (mm)	0.021 (inch)
Cascade L	0.5 – 6 (mm)	37 (mm)	193 (cm)	150 (cm)	4 – 6 (mm)	0.021 (inch)
Cascade L Agile	0.5 – 6 (mm)	21 (mm)	191 (cm)	150 (cm)	4 – 6 (mm)	0.021 (inch)

Table 2 - Braided net's effective length per vessel diameter

Model	Vessel Diameter									
	1.5 mm	2.0 mm	2.5 mm	3.0 mm	3.5 mm	4.0 mm	4.5 mm	5.0 mm	5.5 mm	6.0 mm
Cascade 17	11 mm	9 mm	7 mm	5 mm	---	---	---	---	---	---
Cascade M	---	30 mm	28 mm	25 mm	22 mm	18 mm	---	---	---	---
Cascade M Agile	---	20 mm	17 mm	14 mm	10 mm	7 mm	---	---	---	---
Cascade L	---	---	---	---	---	24 mm	21 mm	18 mm	15 mm	12 mm
Cascade L Agile	---	---	---	---	---	17 mm	14 mm	12 mm	8 mm	6 mm

CAUTION: If any resistance is detected, stop re-sheathing, identify the problem and try again. If the problem persists, use a new device.


Figure 2 - Handle actuation


Figure 3 - Handle modes

INSTRUCTIONS FOR USE

1. Administer anti-coagulation and anti-platelet medications, as per standard practice guidelines.
2. Using angiographic radiography, determine the location and size of the aneurysm to be treated, and the corresponding parent vessel.
3. Select the appropriate Cascade model and size to match the parent vessel and aneurysm properties according to tables 1 and 2.
4. Using conventional catheterization techniques and fluoroscopic visualization, position a suitable guide catheter near the aneurysm's parent vessel. Connect a rotating hemostasis valve (RHV) to the fitting of the guide catheter, and then connect it to the continuous flush line.
5. Connect a second RHV to the fitting of the microcatheter and connect it to the continuous flush line.
6. Using conventional catheterization techniques and fluoroscopic visualization, and aided by a suitable guidewire, advance a microcatheter within the parent vessel and position it at the deployment site.
7. Remove the guidewire from the microcatheter and flush the microcatheter.
8. Closely observe the Cascade and verify that the tip of the device is not protruding beyond the loading tube.
9. Loosen the RHV and introduce the Cascade's loading tube through the port. Seat it firmly in the microcatheter's hub.
10. Tighten the RHV and flush the braided net by allowing 30 seconds of backflow through the loading tube. This will properly flush the braided net and prevent air embolism.
11. Using a smooth and continuous motion, advance the Cascade through the loading tube until the metal portion of the actuation shaft is inserted into the microcatheter. Retract the loading tube from the RHV and retighten if necessary.
12. Advance the Cascade until the fluorosaver markings reach the RHV. Then, aided by fluoroscopic visualization, advance the Cascade until its tip aligns with the tip of the microcatheter.

WARNINGS

- > Do not allow the distal tip to exit the loading tube before fully seating the loading tube within the microcatheter hub.
- > Do not advance the Cascade against resistance without careful assessment of cause. Advancing the device against resistance may result in vessel and/or device damage. If any resistance is met during insertion or when advancing the Cascade, make sure the path through which the Cascade is inserted and advanced is as straight as possible. If resistance is still prominent and cause cannot be determined, withdraw the Cascade completely and try again.
- 13. If needed, reposition the microcatheter in accordance with the vasculature of the aneurysm's parent vessel.
- 14. Unsheathe the Cascade by retracting the microcatheter as needed.
- 15. Using the handle, increase the diameter of the Cascade's braided net to provide mechanical support for the coiling procedure. Make sure the aneurysm neck is covered by the braided net as required.

Note: The braided net can be shortened by partially re-sheathing it with the microcatheter, to comply with unique vasculature anatomy and provide the support required for the specific aneurysm location, position, and shape.

16. Before detaching each coil, carefully observe the braided net to verify that no coil has penetrated the net. If a coil penetrated the net, do not detach the coil. Retract the coil until it is freed from the net and redeploy the coil.

Note: During the embolization procedure, the diameter of the Cascade's braided net can be decreased, and the device can be repositioned, as needed. When doing so, make sure the coils do not herniate into the parent vessel or into the braided net.

17. When the coiling procedure is complete, observe the Cascade's braided net carefully and verify that no coil has penetrated the net.
18. Decrease the diameter of the braided net completely, and re-sheath it with the microcatheter.
19. Retract and remove the Cascade.

WARNINGS

- > Do not apply excessive force when adjusting the cascade diameter. Excessive force may cause oversizing, leading to vessel damage and/or other complications.
- > If the braided net fails to respond to handle manipulations, carefully reposition the Cascade and check for responsiveness. If it continues to fail to respond, carefully withdraw the Cascade for inspection.
- > Do not withdraw the Cascade against resistance without careful assessment of the cause. Withdrawing the Cascade against resistance may result in vessel and/or device damage. If cause cannot be determined, gently re-sheath the cascade with the microcatheter, and remove the device and microcatheter as a single unit.

STORAGE

Store in cool, dry, and dark place.

WARRANTY DISCLAIMER

Although this product has been manufactured under carefully controlled conditions, the manufacturer has no control over the conditions under which this product is used. The manufacturer therefore disclaims all warranties, both expressed and implied with respect to the product, but not limited to, any implied warranty of merchantability or fitness for a particular purpose. The manufacturer shall not be liable to any person or entity for any medical expenses or any direct, incidental or consequential damages caused by any use, defect, failure or malfunction of the product, whether a claim for such damages is based upon warranty, contract, tort or otherwise. No person has any authority to bind the manufacturer to any representation or warranty with respect to the product. The exclusions and limitations set out above are not intended to, and should not be construed so as to contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held illegal, unenforceable or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of this Disclaimer of Warranty shall not be affected and all rights and obligations shall be construed and enforced as if this Disclaimer of Warranty did not contain the particular part or term held to be invalid.

MANUFACTURER

Perflow Medical Ltd.
4 Hatzoran St. Netanya, 4250604 Israel
Tel: +972-3-6544011
Fax: +972-3-5546211
www.perflow.com

AUTHORIZED EUROPEAN REPRESENTATIVE

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MedNet EC-REP GmbH