varian

TANDEM

EMBOZONE TANDEM MICROSPHERES

Microspheres for Embolization

Instructions for Use	2
Instrucciones de uso	20
Mode d'emploi	39
Bedienungsanleitung	58
Istruzioni per l'uso	77

P1064531-001-D 09/2024

TABLE OF CONTENTS

WARNING	3
DEVICE DESCRIPTION.....	3
Syringe	3
User Information.....	3
Contents.....	3
INTENDED USE/INDICATIONS FOR USE.....	3
PERFORMANCE CHARACTERISTICS OF THE DEVICE	4
SIDE EFFECTS AND RESIDUAL RISKS	4
Risk Assessments (Limitations)	4
Residual Risks by Event.....	4
Material Hazard Mitigations.....	4
Chemical Hazard Mitigations	4
Medical-Surgical Hazard Mitigations.....	5
Microspheres Hazard Mitigations.....	5
Device (Normal Function) Hazard Mitigations	5
Labeling and Packaging Hazard Mitigations.....	5
CLINICAL SUMMARY AND BENEFITS	5
SUMMARY OF SAFETY AND CLINICAL PERFORMANCE	5
MAGNETIC RESONANCE IMAGING	5
CONTRAINDICATIONS.....	5
WARNINGS.....	6
PRECAUTIONS.....	6
ADVERSE EVENTS	7
HOW SUPPLIED.....	7
Handling and Storage.....	7
OPERATIONAL INSTRUCTIONS.....	7
Additional Materials and Equipment.....	7
Catheter Compatibility	7
Pre-Procedure Evaluation	8
Procedure Preparation.....	8
Drug Selection and Loading	8
Loading Information.....	8
Additional Materials and Equipment.....	8
Maximum Loading Dose	9
Step 4: Storage (optional)	13
Method 1 - Loading Steps.....	10
Method 2 - Loading Steps.....	11
Method 3 - Loading Steps.....	14
Contrast Mixing and Microsphere Delivery.....	17
POST IMPLANTATION PATIENT INSTRUCTIONS	18
CONTACT VARIAN CUSTOMER SUPPORT/REPORT INCIDENTS	18
WARRANTY	19

TANDEM

EMBOZONE TANDEM MICROSPHERES

Microspheres for Embolization

Rx ONLY

WARNING

Contents supplied STERILE using a steam process. Do not use if sterile barrier is damaged. If damage is found, call your Customer Service representative.

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

DEVICE DESCRIPTION

TANDEM Microspheres are spherical, tightly calibrated, biocompatible, non-resorbable, hydrogel microspheres coated with an inorganic perfluorinated polymer (Polyzene™-F). TANDEM Microspheres are opaque in color, supplied in prefilled syringes, and available in three sizes suitable for embolic therapy. TANDEM Microspheres can be loaded with chemotherapy drugs.

TANDEM Microspheres has a three-layer structure:

- A) Core (98%) Sodium Polymethacrylate (tinted) Barium Sulfate
- B) Core-Shell-Interface (<1%): 3-Aminopropyltrimethoxysilane
- C) Shell (1%): Polyzene®-F

Syringe

TANDEM Microspheres are supplied in a 20 ml syringe prefilled with 2 ml or 3 ml of product suspended in a non-pyrogenic, sterile transport solution of physiological saline. The total volume of TANDEM Microspheres including transport solution is approximately 7 ml. Prefilled syringes of TANDEM Microspheres are packaged in a sterile, sealed tray with a peel-away lid. The label indicates the specific size of the microspheres contained in the syringe.

User Information

Contents

Quantity	Material
1x	TANDEM Microspheres Prefilled Syringe

INTENDED USE/INDICATIONS FOR USE

TANDEM Microspheres are indicated for embolization of blood vessels supplying primary hypervascular tumors or metastases in the liver.

INTENDED USERS AND TARGET PATIENT GROUPS

TANDEM Intended Users:

TANDEM Microspheres are intended to be used by Interventional Radiologists, or other specialized physicians trained in vascular embolization procedures.

Target Patient Groups:

Embolization using TANDEM Microspheres is intended for adult patients with diagnosis suitable for embolotherapy.

PERFORMANCE CHARACTERISTICS OF THE DEVICE

The Embosphere Family Microspheres essential performance objectives are defined as the procedural/technical success rate of embolization of the target vessel (e.g. cessation of or diminished blood flow) and the rates of:

- tumor/lesion response (e.g. overall response rate of complete response (CR) plus partial response (PR)), time to progression and overall survival.

The Embosphere Family Microspheres shall maintain these essential performances:

- ≥95% within the specified size range.
- Compressible up to 30% for smooth delivery through a catheter and robust to facilitate smooth delivery; elasticity allow the microsphere to rebound and maintain its shape.
- No aggregation and minimal fragmentation during delivery.
- Suspended in contrast media and maintains homogenous suspension (post contrast introduction) within the time indicated below:
 - Suspend within 5-10 minutes dependent on size.
 - Stay suspended for more than 1 or 2 minutes dependent on size.
- Suitable for permanent patient implantation.

SIDE EFFECTS AND RESIDUAL RISKS

Risk Assessments (Limitations)

Residual risks are risks that remain after all possible efforts to eliminate them. A Low risk is broadly acceptable (for example, results in temporary discomfort). A Medium risk is acceptable only if the benefit of the device outweighs the residual risk, and the residual risk is reduced as low as reasonably achievable.

No identified single fault or use errors remain that could result in death or serious injury when a clinician of average training and experience uses the product as directed. Clinician of average training and experience must employ industry-standard or higher, or consensus quality assurance and safety measures appropriate for the treatments planned and delivered.

Residual Risks by Event

The lowest attainable risk varies by hazardous event.

TABLE A - Residual Risks by Event

HAZARD	LOWEST ATTAINABLE RISK
Materials Infectious substances, biological products or biologically incompatible materials	Medium
Chemical Overdosing, Underdosing, Drug Release	Low
Medical-Surgical Pathogens or Other Toxicants	Low
Microspheres Size, Delivery, Debris	Low
Device Normal Function	Low
Labeling and Packaging	Low

Material Hazard Mitigations

Varian has defined cleaning procedures and conducted biocompatibility assessments for all parts that come into contact with patients and clinical personnel. Clinicians should adhere to procedures and follow good clinical practice set by their clinic or the hospital environment.

Chemical Hazard Mitigations

Clinicians may control risks by procedures and following good clinical practice set by their clinic or the hospital environment; this includes following the Instructions for Use.

Medical-Surgical Hazard Mitigations

Clinicians may control risks by procedures and following good clinical practice set by their clinic or the hospital environment; this includes following the Instructions for Use.

Microspheres Hazard Mitigations

Varian has defined and conducted assessments on Microspheres that comply to industry standards, that mitigate any further risks related to device functionality.

Device (Normal Function) Hazard Mitigations

Clinicians may control risks by procedures and following good clinical practice set by their clinic or the hospital environment; this includes following the Instructions for Use.

Additionally, Varian has defined and conducted assessments on Microspheres that comply to industry standards, thus mitigating any further risks related to device functionality.

Labeling and Packaging Hazard Mitigations

Varian has defined and conducted assessments on Microsphere's packaging and labeling that comply to industry standards, thus mitigating any further risks related to device functionality.

CLINICAL SUMMARY AND BENEFITS

The intended clinical benefit is occlusion of blood vessels which leads to reduction or occlusion of blood flow ultimately leading to necrosis of the target tissue. The goal is the reduction or cessation of growth of the target tissue or reduction of the patient's symptoms. Typical target tissues are benign or malignant tumors, arteriovenous malformations or bleeding vessels other than the central nervous system. Untreated target tissues would typically continue to grow or cause worsening symptoms and ultimately have a negative impact on a patient's health, including death.

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE

The Summary of Safety and Clinical Performance (SSCP) is intended to provide public access to an updated summary of the main aspects of the safety and clinical performance of this device. The SSCP is not intended to replace the Instructions for Use as the main document to ensure the safe use of the device, nor is it intended to provide diagnostic or therapeutic suggestions to intended users or patients.

The SSCP can be found at www.varian.com/SSCP/TANDEM

MAGNETIC RESONANCE IMAGING

TANDEM Microspheres are MR Safe.

CONTRAINDICATIONS

Embolization procedures shall not be performed if:

- Patient is unable to tolerate vascular occlusion procedures.
- Vascular anatomy precludes correct catheter placement or embolic injection.
- Presence or likely onset of vasospasm.
- Presence of a blood coagulation disorder that would prohibit arterial punctures.
- Presence of severe atheromatous disease that would preclude correct catheter placement.
- Presence of patent extra-to-intra-cranial anastomoses or shunts from the arterial to the venous circulation.
- Presence of collateral vessel pathways which could potentially endanger non-targeted tissue during an embolization procedure.
- Presence of any vasculature where TANDEM Microspheres could pass directly into the central nervous system, central circulatory system or other non-target territories.
- Patient has high-flow arteriovenous shunt with diameter greater than the selected TANDEM Microspheres.
- Patient is pregnant.
- Patient has known allergies to barium sulfate, 3-aminopropyltrialkoxysilane, polyphosphazene, IV radiopaque contrast agent.

WARNINGS

Vascular embolization is a high risk procedure. The procedure should be performed by specialized physicians experienced in vascular embolization and chemoembolization procedures:

- Risks of radiation from angiography and fluoroscopy used to visualize the blood vessels during embolization, which may include a radiation burn and risks to future fertility.
 - Do not use TANDEM Microspheres in conjunction with embolization devices based on organic solvents such as ethyl alcohol or dimethyl sulfoxide (DMSO) at the same embolization site.
 - Do not use ionic contrast agent with this product. Ionic contrast agents could alter the microsphere characteristics resulting in microsphere deformation and procedure failure.
 - Do not use heparinized saline as this could lead to microsphere agglomeration. Agglomeration may impede microsphere delivery through the catheter or result in non-target embolization.
 - Should catheter obstruction occur, remove the catheter from the patient. Do not use forceful injection, guidewires or other instruments to dislodge the blockage.
-

PRECAUTIONS

To maintain safety, the following precautions shall be considered:

- Each package of TANDEM Microspheres is intended for single patient use only. Discard any unused material.
- Physicians using TANDEM Microspheres should have appropriate training and experience in a related interventional procedure.
- Do not use TANDEM Microspheres if the sterile barrier, the syringe or the package appears to be opened or damaged prior to use.
- Do not use TANDEM Microspheres that have been improperly stored or mishandled.
- The physician should carefully select the size and quantity of TANDEM Microspheres according to the lesion to be treated based on the physician's experience and currently available scientific evidence.
- Physicians must decide the most appropriate time to stop the infusion of TANDEM Microspheres. Typically the artery will accept fewer TANDEM Microspheres as the treatment progresses. Proximal slowing or termination of flow may indicate that the vessel or the target area is occluded by TANDEM Microspheres. Careful fluoroscopic monitoring is required.
- Microsphere embolization must be performed slowly. The injection speed and manner must be controlled. Excessive injection rate may result in retrograde flow in the vessel leading to embolization of other non-target healthy tissue or organs.
- If arteriovenous anastomoses, branch vessels which lead away from the targeted embolization area, or emergent vessels not evident prior to embolization are present, it can lead to non-target embolization and cause severe complications for the patient.
- Ischemia of tissue adjacent to the targeted area may result from post-embolization swelling. Therefore, special care should be taken to avoid such ischemia of nontolerant, non-targeted tissue such as the nervous system.
- Consider upsizing TANDEM Microspheres if angiographic appearance of embolization does not quickly appear during injection of the microspheres.
- If there are any symptoms of non-target embolization during injection, consider stopping the procedure to evaluate the possibility of shunting. Such symptoms may include changes in patient vital signs, such as hypoxia or central nervous system changes.

Precautions specific to use of small microspheres:

- TANDEM Microspheres smaller than 100 μm can migrate to distal anastomotic feeders and embolize circulation to distal tissue. For this reason, smaller particles have a greater likelihood of causing unwanted ischemic injury. This should be considered prior to starting the embolization procedure.

Interaction with Pharmaceuticals:

There are no known chemical interactions between TANDEM Microspheres and pharmaceuticals.

ADVERSE EVENTS

Potential adverse events associated with the use of TANDEM Microspheres include, but may not be limited to:

- Allergic reaction
- Capillary bed saturation and tissue damage
- Cerebrovascular accident (CVA)
- Complications related to catheterization (e.g., hematoma at the site of entry, clot formation at the tip of the microcatheter and subsequent dislodgement, vasospasm, nerve injury, vessel trauma [e.g., dissection, perforation, rupture])
- Death
- Foreign body reactions (e.g., pain, rash, fever, inflammation)
- Hemorrhage
- Incomplete occlusion of vascular beds or territories may give rise to the possibility of post-procedural hemorrhage, development of alternative vascular pathways, recanalization, or recurrence of symptoms
- Infection
- Ischemia at an undesirable location
- Ischemic infarction
- Neurological deficits
- Post-embolization syndrome
- Pulmonary embolization
- Thrombosis
- Undesirable reflux, passage/migration or placement of TANDEM Microspheres, resulting in non-target embolization
- Vessel or lesion rupture

HOW SUPPLIED

Prefilled syringes of TANDEM Microspheres are packaged in a sealed tray with a peel-away lid.

Sterile

Non-pyrogenic

Do not use if package is opened or damaged.

Do not use if labeling is incomplete or illegible.

Handling and Storage

Store in a cool, dry, dark place.

Product must be used prior to expiration date on label.

OPERATIONAL INSTRUCTIONS

Additional Materials and Equipment

(1x) Microsphere injection syringe

(1x) 3-way stopcock

Delivery catheter

Non-ionic contrast media

Water for Injection (WFI)

Catheter Compatibility

TANDEM Microspheres are designed to be used with a variety of catheters and microcatheters. Select a delivery catheter of appropriate size, suitable for the dimensions of the target vessels. TANDEM Microspheres can tolerate temporary compression to facilitate passage through the delivery catheter. Utilize the catheter's minimum inner diameter measurement to determine catheter-to-microsphere compatibility. For 75-100 μ m Microsphere sizes, the compatible catheter size minimum ID is 0.008 inches (0.2mm).

Pre-Procedure Evaluation

Prior to the embolization procedure, evaluate the vascular anatomy associated with the lesion using high resolution imaging.

Procedure Preparation

1. Select the size of TANDEM Microspheres according to the size of the vessel identified.

Caution: TANDEM Microspheres smaller than 100 μm can migrate to distal anastomotic feeders and embolize circulation to distal tissue. For this reason, smaller particles have a greater likelihood of causing unwanted ischemic injury. This should be considered prior to starting the embolization procedure.

2. Verify that the sterile packaging was not compromised.

Drug Selection and Loading

TANDEM Microspheres can be loaded with chemotherapeutic drugs. When used for drug loading, drug loading should be done under a physician's direction, choice, and responsibility.

The Loading Information that follows will provide guidelines for the loading of:

- Doxorubicin hydrochloride (lyophilized powder)
- Irinotecan hydrochloride (20 mg/ml solution)
- Doxorubicin hydrochloride (2 mg/ml solution)

The loading of TANDEM Microspheres with doxorubicin hydrochloride (lyophilized powder) containing the following has been evaluated:

- Non-ionizable excipient(s): lactose monohydrate and/or methyl-4-hydroxybenzoate

The loading of TANDEM Microspheres with irinotecan hydrochloride (20 mg/ml solution) containing the following has been evaluated:

- Non-ionizable excipient(s): sorbitol powder
- pH adjuster(s): sodium hydroxide, hydrochloric acid, lactic acid
- Water for injection

The loading of TANDEM Microspheres with doxorubicin hydrochloride (2 mg/ml solution) containing the following has been evaluated:

- Excipient(s): sodium chloride (up to 0.9% or 9 mg/ml)
- pH adjuster(s): sodium hydroxide and/or hydrochloric acid
- Water for injection

The loading of TANDEM Microspheres with other drugs has not been evaluated.

Caution: Do not use liposome formulations of the drugs as it may interfere with drug-loading capability.

Loading Information

Refer to the drug manufacturer's package insert for drug specific handling, contraindications, warnings, precautions, and adverse events. In order to minimize the risk of microbiological contamination, TANDEM Microspheres' drug-loading should be performed under controlled aseptic conditions.

Additional Materials and Equipment

Water for injection (WFI)

(2x) 5 μm sterile filter needle

18 gauge - 20 gauge sterile needle

Waste container (e.g. empty 100 ml solution bag)

Delivery syringe (20 ml syringe)

Loading syringe (10 ml, 20 ml, 60 ml or 100 ml; see Maximum Loading Dose for the appropriate size)

Doxorubicin hydrochloride (lyophilized powder), Irinotecan hydrochloride (20 mg/ml solution), or Doxorubicin hydrochloride (2 mg/ml solution)

Maximum Loading Dose

The maximum loading dose is 50 mg doxorubicin/1 ml microspheres or 50 mg irinotecan/1 ml microspheres and can be achieved using the following guidelines.

TABLE B**Summary of drug selection and loading methods.**

DRUG	MICROSPHERE SYRINGE SIZE	TOTAL VOLUME OF DRUG SOLUTION	DOSE	LOADING SYRINGE	LOADING STEPS PAGE #
METHOD 1					
Doxorubicin Powder, 20 mg/ml	2 ml	Reconstitute with 5 ml of WFI	100 mg	≥ 10 ml	10
Doxorubicin Powder, 20 mg/ml	3 ml	Reconstitute with 7.5 ml of WFI	150 mg	≥ 20 ml	10
Irinotecan, 20 mg/ml	2 ml	5 ml	100 mg	≥ 10 ml	10
Irinotecan, 20 mg/ml	3 ml	7.5 ml	150 mg	≥ 20 ml	10
METHOD 2					
Doxorubicin Solution, 2 mg/ml	2 ml	50 ml	100 mg	≥ 60 ml	11
Doxorubicin Solution, 2 mg/ml	3 ml	75 ml	150 mg	100 ml	11
METHOD 3					
Doxorubicin Solution, 2 mg/ml	3 ml	75 ml	150 mg	2x 60 ml	14

TABLE C**Minimum loading times with achievable drug loading percentages.**

DESCRIPTION	MINIMUM LOADING TIME	LOADING PERCENTAGE
40µm Powdered Doxorubicin	60 minutes	75% - 100%
75µm Powdered Doxorubicin	60 minutes	96% - 100%
100µm Powdered Doxorubicin	60 minutes	96% - 100%
40µm Liquid Doxorubicin	120 minutes	75% - 100%
75µm Liquid Doxorubicin	120 minutes	96% - 100%
100µm Liquid Doxorubicin	120 minutes	96% - 100%
40µm Irinotecan	30 minutes	96% - 100%
75µm Irinotecan	30 minutes	96% - 100%
100µm Irinotecan	30 minutes	96% - 100%

Caution: Do not exceed maximum loading dose.

METHOD 1 - LOADING STEPS

Doxorubicin Powder & Irinotecan Solution

- Doxorubicin 20 mg/ml solution reconstituted from powder
- Irinotecan 20 mg/ml solution

Step 1: Microspheres Preparation

- 1.1 Open box and remove inner packaging. Verify sterile package is free of damage. Do not use if sterile barrier is damaged.
- 1.2 Peel lid to open sterile package and remove syringe from tray.
- 1.3 Gently swirl contents and allow microspheres to settle towards the plunger by standing syringe upright. It may take several minutes to fully settle.
- 1.4 While keeping syringe upright and taking care not to disrupt settled microspheres, attach a sterile filter needle to the syringe. With the filter needle pointed upward, eject transport solution (Figure 1). Ensure enough transport solution remains inside the syringe to cover the microspheres.
- 1.5 Remove the filter needle and attach a sterile needle.


Figure 1.

Step 2: Drug-Loading

- 2.1 Transfer required amount of drug from its container to the syringe by drawing solution into syringe via the attached sterile needle (Figure 2a). Leave a small air bubble (1 ml or less) in the syringe for better mixing.


Figure 2a.

- 2.2 Re-cap syringe.
- 2.3 Mix by gently inverting/swirling the syringe every 5 minutes (Figure 2b). Hard shaking may cause foam.


Figure 2b.

2.4 Allow for sufficient loading time. Recommended minimum loading time to achieve drug-loading:

- Doxorubicin (solution, 20 mg/ml, reconstituted from powder) minimum loading time is 60 minutes.
- Irinotecan (solution, 20 mg/ml) minimum loading time is 30 minutes.

Step 3: Storage (optional)

After drug-loading is complete, if not used immediately, store in the dark at 2°C - 8°C and use within 48 hours. Please refer to the drug manufacturer's package insert for specific drug information.

Step 4: Supernatant Removal

- 4.1 Prior to use, allow the drug-loaded microspheres to settle towards the plunger by standing syringe upright. It may take several minutes to fully settle.
- 4.2 While keeping syringe upright and taking care not to disrupt settled microspheres, attach a sterile filter needle to the syringe. With the filter needle pointed upward, eject supernatant (Figure 3). Ensure enough supernatant remains inside the syringe to cover the microspheres.


Figure 3.

4.3 Re-cap syringe.

4.4 Go to page 17: Contrast Mixing and Microsphere Delivery.

METHOD 2 - LOADING STEPS

Doxorubicin Solution Using One Syringe

- Doxorubicin 2 mg/ml solution

- For 2 ml of microspheres, use 50 ml (100 mg) of drug solution in a 60 ml syringe.

- For 3 ml of microspheres, use 75 ml (150 mg) of drug solution in a 100 ml syringe. If a 100 ml syringe is not available see page 14 for loading instructions.

Step 1: Microspheres Preparation (including Transfer into Loading Syringe)

- 1.1 Open box and remove inner packaging. Verify sterile package is free of damage. Do not use if sterile barrier is damaged.
- 1.2 Peel lid to open sterile package and remove product syringe from tray.
- 1.3 Select an appropriately sized loading syringe (60 ml or 100 ml). Connect product syringe and loading syringe to a 3-way stopcock (Figure 4).
- 1.4 Gently swirl the product syringe to suspend microspheres in transport solution. Transfer microspheres from product syringe to loading syringe via attached 3-way stopcock.
- 1.5 Disconnect loading syringe from 3-way stopcock.
- 1.6 Gently swirl contents and allow microspheres to settle towards the plunger by standing syringe upright. It may take several minutes to fully settle.
- 1.7 While keeping syringe upright and taking care not to disrupt settled microspheres, attach a sterile filter needle to the loading syringe. With the filter needle pointed upward, eject transport solution (Figure 4). Ensure enough transport solution remains inside the loading syringe to cover the microspheres.
- 1.8 Remove the filter needle and attach a sterile needle.


Figure 4.

Step 2: Drug-Loading

- 2.1 Transfer required amount of drug from its container to the loading syringe by drawing solution into loading syringe via the attached sterile needle (Figure 5a). Leave a small air bubble (1 ml or less) in the loading syringe for better mixing.
- 2.2 Re-cap loading syringe.
- 2.3 Mix by gently inverting/swirling the loading syringe every 5 minutes (Figure 5b). Hard shaking may cause foam.
- 2.4 Allow for sufficient loading time. Recommended minimum loading time to achieve drug-loading for doxorubicin (solution, 2 mg/ml) is 120 minutes.


Figure 5a.


Figure 5b.

Step 3: Transfer into 20 ml delivery syringe

- 3.1 Gently swirl contents and allow microspheres to settle towards the plunger. It may take several minutes to fully settle.
- 3.2 Attach a sterile filter needle to the loading syringe. With the filter needle pointed upward, begin to eject transport solution, leaving ~15 ml of the transport solution in the loading syringe.
- 3.3 Connect loading syringe and delivery syringe to a 3-way stopcock.
- 3.4 Gently swirl to suspend the microspheres in the drug solution in the loading syringe. Transfer microspheres from loading syringe to delivery syringe via attached 3-way stopcock.

Step 4: Storage (optional)

After drug-loading is complete, if not used immediately, store in the dark at 2°C - 8°C and use within 48 hours. Please refer to the drug manufacturer's package insert for specific drug information.

Step 5: Supernatant Removal

- 5.1 Prior to use, allow the drug loaded microspheres to settle towards the plunger by standing syringe upright. It may take several minutes to fully settle.
- 5.2 While keeping syringe upright and taking care not to disrupt settled microspheres, attach a sterile filter needle to the delivery syringe. With the filter needle pointed upward, eject supernatant (Figure 6). Ensure enough supernatant remains inside the delivery syringe to cover the microspheres.


Figure 6.

5.3 Re-cap delivery syringe.

5.4 Go to page 17: Contrast Mixing and Microsphere Delivery

METHOD 3 - LOADING STEPS

Doxorubicin Solution Using Two Syringes

- Doxorubicin 2 mg/ml solution

- For 3 ml of microspheres, use 75 ml of drug solution (150mg) in 2 x 60 ml syringes

Step 1: Microspheres Preparation

- 1.1 Open box and remove inner packaging. Verify sterile package is free of damage. Do not use if sterile barrier is damaged.
- 1.2 Peel lid to open sterile package and remove product syringe from tray.
- 1.3 Gently swirl and allow microspheres to settle towards the plunger by standing syringe upright. It may take several minutes to fully settle.
- 1.4 While keeping syringe upright and taking care not to disrupt settled microspheres, attach a sterile filter needle to the product syringe. Insert the filter needle into a waste container, with the needle pointed upward, eject transport solution until a small amount (1 ml or less) of solution remains in the product syringe. Do not compress the microspheres by leaving a thin layer of transport solution covering the microspheres. (Figure 7)


Figure 7.

- 1.5 Remove the filter needle and recap the product syringe.

Step 2: Drug-Loading

- 2.1 Draw up the drug using 2 x 60 ml loading syringes, 37.5 ml of drug in each syringe (75 mg total). Leave a small air bubble (1 ml or less) in each of the loading syringes for better mixing. (Figure 8).


Figure 8.

- 2.2 Connect product syringe and 1 of the 60 ml loading syringes to a 3-way stopcock.
- 2.3 Gently push some of the drug solution into the microsphere syringe and then continue to slowly push the microsphere-drug solution back and forth a few times between both syringes so that the microspheres are adequately suspended. Push the entire microsphere-drug solution into the 60 ml syringe (Figure 9).


Figure 9.

- 2.4 Attach the microsphere-drug loading syringe and the second drug syringe to a 3-way stopcock (Figure 10).


Figure 10.

- 2.5 Gently mix by slowly pushing the contents from one syringe to the other. Mix solution until microspheres are suspended (Figure 11).


Figure 11.

- 2.6 Adjust volume of syringes into equal amounts. Each syringe should have a total volume of 39 ml (37.5 ml drug solution + 1.5 ml of microspheres). Leave a small air bubble (1 ml or less) in each syringe for better mixing.
- 2.7 Re-cap each of the two loading syringes.
- 2.8 Mix by gently inverting/swirling the loading syringe every 5 minutes (Figure 12). Hard shaking may cause foam.
- 2.9 Allow for sufficient loading time. Recommended minimum loading time to achieve drug-loading for doxorubicin (solution, 2 mg/ml) is 120 minutes.


Figure 12.

Step 3: Storage (optional)

- 3.1 After drug-loading is complete, if not used immediately, store the syringes in the dark at 2°C - 8°C and use within 48 hours. Please refer to the drug manufacturer's package insert for specific drug information.

Step 4: Supernatant Removal

- 4.1 Gently swirl contents and allow microspheres to settle towards the plunger. It may take several minutes to fully settle (Figure 13).


Figure 13.

- 4.2 Attach a sterile filter needle to one of the loading syringes. With the filter needle pointed upward, expel supernatant allowing ~10 ml to remain in the syringe. Repeat step with the 2nd loading syringe (Figure 14).


Figure 14.

- 4.3 Recap the syringes.

Step 5: Transfer into 20ml delivery syringe

- 5.1 Connect one loading syringe and 20 ml delivery syringe to a 3-way stopcock.

- 5.2 Transfer microspheres from one loading syringe to the delivery syringe via attached 3-way stopcock. Disconnect the 1st loading syringe and attach the 2nd loading syringe. Transfer the microspheres to the delivery syringe (Figure 15).


Figure 15.

- 5.3 Gently swirl and allow microspheres to settle towards the plunger in the 20 ml delivery syringe.
- 5.4 Attach filter needle to the 20 ml syringe.
- 5.5 Insert needle into a waste container and expel excess supernatant until a small amount (>1 ml) of supernatant remains in the 20 ml delivery syringe. Do not compress the microspheres by leaving a thin layer of supernatant covering the microspheres in the syringe (Figure 16).


Figure 16.

- 5.6 Recap the syringe.

CONTRAST MIXING AND MICROSPHERE DELIVERY

1. Position the catheter at the desired site and perform baseline angiography to evaluate the blood supply to the lesion.
2. Prior to adding contrast agent to the drug-loaded microspheres, verify that the supernatant has been removed. If not, refer to Supernatant Removal, above.
3. Use only non-ionic contrast agent in accordance with the contrast agent labeling with respect to dosage.

Warning: Do not use ionic contrast agent with this product. Ionic contrast agents could alter the microsphere characteristics resulting in microsphere deformation and procedure failure.

4. Add an appropriate amount of contrast agent to the 20 ml syringe containing microspheres to obtain homogeneous suspension and fluoroscopic visibility. Pure contrast agent or a mixture of contrast agent and WFI may be used.

Caution: Do not use saline as it may interfere with drug-loading capability.

A mixture is especially recommended when less than 50 mg of drug per ml of microspheres has been loaded or when using contrast agent with iodine concentrations higher than 300 mg iodine/ml. The suspension is usually obtained with a mixture of 50% contrast agent and 50% WFI. Contrast agent and WFI can be added in the same proportions to obtain a more dilute suspension. Agitate the microspheres and contrast agent until a homogeneous suspension is achieved.

5. Confirm suspension prior to delivery. If microspheres have begun to settle, agitate to re-suspend before delivery.
6. Purge all air from the 20 ml syringe containing microspheres.
7. Attach the 20 ml syringe containing microspheres to one port of the 3-way stopcock and an injection syringe (1 ml - 3 ml) to another port of the 3-way stopcock. Attach a delivery catheter to the remaining port on the 3-way stopcock. Ensure the 3-way stopcock is securely attached.
8. Draw the microsphere mixture slowly and gently into the injection syringe to minimize the potential of introducing air into the system.
9. Under continuous fluoroscopic control, slowly infuse microspheres into the blood stream. Always inject under free flow conditions by fully opening the 3-way stopcock. To optimize injection through the catheter, it is recommended that the injection syringe remains in a horizontal position during injection.

Caution: Microsphere embolization must be performed slowly. The injection speed and manner must be controlled. Excessive injection rate may result in retrograde flow in the vessel leading to embolization of non-target healthy tissue or organs.

10. Continue infusion until the desired devascularization is achieved.
11. Once the procedural endpoint is reached, wait for 5 minutes to observe whether the microspheres redistribute themselves and re-establish flow to the target. If flow is re-established, inject an additional volume of microspheres until the final procedural endpoint is achieved.
12. At the end of the infusion, remove the catheter while maintaining gentle aspiration to avoid dislodging any residual microspheres still inside the catheter.
13. Discard any opened unused syringes of TANDEM Microspheres.

POST IMPLANTATION PATIENT INSTRUCTIONS

After Microsphere administration, the physician should give the patient the implant card with the appropriate information filled in. This information includes patient name, date of implant, physician name, device name, device model number, batch/lot/serial number. Information as part of the implant card contains legal manufacturer name and address, and MR compatibility, and latex free. Physician should advise the patient to keep the card in a safe place, e.g. his or her wallet, for future reference. The physician should advise the patient that this implant card contains important information related to the Embozene Microsphere Family implant and that the card should be shown to their current and future health care providers.

In the instance of a patient losing their card, the patient should be told to contact Varian Support.

CONTACT VARIAN CUSTOMER SUPPORT/REPORT INCIDENTS

For customer support, contact us at vis.support@varian.com or call +1 (877) 457-2009 or contact your local representative.

Any serious incident that has occurred while using the device should be reported to the manufacturer and the competent authority of the Member State in which the user or patient is established.

To report incidents, contact us at microspherecomplaints@varian.com or call +1 (877) 457-2009 or contact your local representative.

WARRANTY

The Manufacturer warrants that reasonable care has been used in the design and manufacture of this instrument. **This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose.** Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond the Manufacturer's control directly affect the instrument and the results obtained from its use. The Manufacturer's obligation under this warranty is limited to the repair or replacement of this instrument and the Manufacturer shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this instrument. The Manufacturer neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. **The Manufacturer assumes no liability with respect to instruments reused, reprocessed or resterilized and makes no warranties, express or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instruments.**

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All rights reserved.
Varian Medical Systems, Inc
Siemens Healthineers Intl. AG
© 2025 Varian Medical Systems, Inc. or its affiliates.