
Mitral Valve Solutions

Chord-X[®] Mitral Chordal Replacement Products


On-X[®] Mitral Heart Valves





Pre-Measured Chord-X[®] Loops


For Mitral Valve Repair



Chord-X is an innovative tool for mitral valve repair, and is designed to:

- Save time*
- Allow reproducible results
- Simplify mitral valve repair procedures *Compared to artificial chords made by surgeon.

Procedural Overview











Attachment

Implant

Complete

Mitral Valve Repair

- Mitral valve repair is strongly preferred over replacement for primary MR whenever anatomically feasible.¹
- Short- and long-term outcomes of successful valve repair for primary MR exceed those for valve replacement across all age ranges.²
- Chordal replacement may be associated with greater freedom from reoperation and may lead to improved postoperative left ventricular function compared with leaflet resection.³
- Chordal implantation keeps the door open for further treatments in case of MR recurrence.⁴




Chord-X Pre-Measured Loops provide an innovative solution to simplify and standardize mitral valve repair.

If mitral valve repair is not possible, consider your options for mitral valve replacement

Mitral Valve Replacement

Advantages of Mechanical Valves:

- "...mechanical valve prostheses remain the gold standard for valve replacement in younger individuals" when mitral valve repair is not possible.⁵
- Mechanical mitral valve offers survival benefit at 15 years for patients <70⁶ (see graph below)



Disadvantages of Tissue Valves:

- Calcification of leaflets increase valve gradients until reoperation is required⁶
- No guarantee patient will not need long term anticoagulation⁷
- 3-4 times more likely to require reoperation when compared to mechanical valves⁶
 - o Note: 30-day mortality for all reoperative mitral valve replacement is 14%⁶

The On-X Mitral Valve offers the "gold standard" solution for mitral valve replacement.⁵

See Instructions for Use for full details.²




For Mitral Valve Replacement



- Pure pyrolytic carbon for reduced thrombogenicity⁸
- < 4.4 mm Hg across all valve sizes 25 – 33 mm⁹
- Sewing cuff designs made from polytetrafluoroethylene (PTFE)⁹
- Leaflet guards designed to prevent encroachment of subvalvular apparatus

90° Leaflet Opening



Ordering Information

Chord-X Pre-Measured Loops and Adjustable Suture System

ePTFE Suture Diameter: USP 2-0 (Diameter similar to GORE-TEX [®] CV-4)	18mm Taper Point Needle Options, Catalog Numbers	
Loop Length	3/8 Circle	1/2 Circle
Adjustable	CXL-20-1838-0	CXL-20-1812-0
12mm	CXL-20-1838-12	CXL-20-1812-12
16mm	CXL-20-1838-16	CXL-20-1812-16
20mm	CXL-20-1838-20	CXL-20-1812-20
24mm	CXL-20-1838-24	CXL-20-1812-24

1 box = 5 sterile Chord-X packaged units.

Each unit contains a Chord-X prosthesis featuring one suture pair for the papillary muscle and three independent
 pre-measured chord-X prosthesis featuring one suture pair for the papillary muscle and three independent
 pre-measured chord-X prosthesis featuring one suture strands for leaflet attachment.
 Each prosthesis is positioned on disposable foam base-plates for suture management.
 Each unit has double-armed taper point needles, two PTE pledgets: one attached and one independent.

Chord-X Chordal Sizer

Package Quantit		Use	Catalog Number
1 box = 5 indiv packaged st Chord-X Ch Sizers	erile	Single use, disposable device	CXCS
Leaf Hook	►) F	Rotatable 360°	•••• 0. V (B)
Papillary Hook	(Fig 2)	Loop Length Adjustable (Fi 12mm 16mm 20mm 24mm	g 1)

On-X Mitral Valve Standard **Sewing Cuff**



On-X Mitral Valve Conform-X **Sewing Cuff**



Valve Size	Product Code
23*	ONXM-23
25	ONXM-25
27/29	ONXM-27/29
31/33	ONXM-31/33

Valve Size	Product Code
25/33	ONXMC-25/33

Learn more at: www.CryoLife.com

- Nishimura RA, Otto CM, Bonow RO, et al. 2014 AHA/ACC guideline for the management of patients with valvular heart disease: executive summary: a report of the American College of Cardiology/American Heart Association Task Force on Practice Guidelines. J Am Coll Cardiol. 2014;63:2438–88.
- Vasileva CM, Mishkel G, McNeely C, et al. Long term survival of patients undergoing mitral valve repair and replacement: a longitudinal analysis of Medicare fee-forservice beneficiaries. Circulation. 2013;127: 1807–6.
 J Thorac Cardiovasc Surg. 2018 Jan;155(1):120-128.e10. doi: 10.1016/j.jtoxs.2017.07.078. Epub 2017 Aug 24. Systematic review and meta-anal-ysis of chordal replacement versus leaflet resection for posterior mitral leaflet prolapse.
- Mindr V. Desia, Francesco Gragioni, Marco Di Eusnio, Matteo Saccocci, Mauricio Taramasso, Francesco Malsanod, Rakesh M. Suri, A. Marc Gillinov. Outcomes in Degenerative Mitral Regurgitation: Current State-of-the Art and Future Directions. Progress In Cardiovascular Diseases 60 (2017) 370-385

5. Chiu P, Goldstone AB, Woo YJ. Eur J Cardiothorac Surg 2018;54:622-6

6. Goldstone AB et al. N Engl J Med 2017;377:1847-1857

7. Nishimura et al., Circulation. 2017;135:e1159-e1195.

LaGrange Let al., Compatibility of carbon and blood. Hegyeli RJ, Editor, Proceedings, Artificial Heart Program Conference, Washington, DC, June 9-13, 1969, 47-58.

9. On-X Prosthetic Heart Valve Instructions for Use.

9. On-X Horstheic Heart Valve Instructions for Use.
On-X Life Technologis, Inc. is a wholly owned subsidiary of CryoLife, Inc. The On-X Prosthetic Heart Valves and Mitral Chordal Replacement Products are manufactured by On-X Life Technologies, Inc. CryoLife, the snowlake design, Life Restoring Technologies, On-X and Chord-X are registread trademarks owned by CryoLife, Inc. Or its subsidiares. All other registread trademarks are evened by their respective owners. Or 2018 CryoLife, Inc. All rights reserved. Surgical technique is at the discretion of the surgeon, Variations in technique and practices will inevitably and appropriately occur when clinicians take into account the needs of the individual patients, available resources, and limitations unique to an institution or type of practice.

