

Spotlight OPS* Procedure Kits

Spotlight OPS Procedure Kits eliminate the need for a separate dilator and provide a cost-effective alternative to other products.

Benefits

- Tungsten-loaded sheath tip increases echogenicity, allowing for bright and easy tip placement, while the smooth, soft, atraumatic tip and low profile make it easy to navigate.
- Unique self-dilation technology eliminates the need for a separate dilator, resulting in fewer components, fewer exchanges, and potentially faster procedures.

SPOTLIGHT OPS PROCEDURE KITS

Description	Part #	Kit Length
Procedure Kit with Spotlight OPS	H787EVL25OPS5	25 cm
Procedure Kit with Spotlight OPS	H787EVL55OPS5	55 cm
Procedure Kit with Spotlight OPS	H787EVL80OPS5	80 cm



VenaCure^{Ops}
Spotlight One-Piece Sheath Fiber